SmartPort⁺

A Patient's Guide

	English.....	1
	Español.....	15
	Ελληνικά	31
	Polski	47
	Italiano	63

IMPORTANT NOTES:

- The safety measures described in the patient literature, product labeling, and instructions for use are designed to minimize risks of port-related complications.
- Any serious incident which has occurred with the use of this device should be reported to your healthcare provider immediately. It is also important for product-related concerns or complaints to be shared with AngioDynamics at complaints@angiodynamics.com
- Product-related concerns or complaints can be submitted by patients and/or healthcare providers.
- This guide is meant to be educational. It is not meant to replace medical care or a doctor's advice. Please talk to your healthcare provider if you have questions or concerns.
- For a copy of the Summary of Safety and Clinical Performance for these devices, please review Eudamed (ec.europa.eu/tools/eudamed) referencing Basic UDI-DI # 50516840024U or contact AngioDynamics Customer Service at 1-800-772-6446.
- More information about your SmartPort⁺ port is available by calling the AngioDynamics Vascular Access Information Line 1-800.772.6446

Important instructions about Implant Card

- You will receive an Implant Card with your port.
- Always carry the Implant Card with you.
- The Implant Card has important information for your healthcare provider.

angiodynamics

www.angiodynamics.com

AngioDynamics, Inc.
 603 Queensbury Avenue
 Queensbury, NY 12804 USA
 USA Customer Service 800-772-6446

Implantable Port System / Sistemas de puertos implantables / Systèmes de chambres implantables / Implantierbare Portsysteme / Sistemi di porte impiantabili / Implanterbare poortsystemen / Implanterbare portsystemer / Implanterbara portsystem / Συστήματα εμφυτεύσιμης θύρας / Sistemas de Porta Implantável / Beültethető portrendszer / Systémy implantabilních portů / Wszczepialne systemy portów / Implanterbare port-systemer / 이식형 포트 시스템 / Implant Edilebilir Port Sistemleri / Implantuojamos prievadų sistemos / Implanteritavad pordisüsteemid / Systeme de porturi implantabile / Systémy implantovatelných portov / Implantabilni port sustavi / Системы на имплантируемые порты / Sistemi implantabilnih portova / Имплантируемые инфузионные системы / Implanitovava portijärjestelmä / أنظمة المنفذ القابلة للزرع / Implantējamās porta sistēmas / Имплантируемые системы от портов / Sistemi vsadne Valvule

SmartPort⁺

H787CT50LTBANFRVIO

UDI-DI: (01)15051684017473

LOT

55515551

T1990000201 A

Table of Contents

Introduction.....	4
What are the SmartPort ⁺ implantable ports made of?	4
How is the port placed?	5
What to expect during recovery from port placement	6
How is the port accessed and used?.....	6
How to identify the SmartPort ⁺ device	6
What is a “CT” or “CECT” Scan, and how is the port used during these types of tests	7
How long can a port remain in place?	7
Commonly asked questions	8
Important information to report to your healthcare provider	9
Patient checklist.....	9
Important Notes and Instructions from the Healthcare Provider:	10
Commonly used terms.....	13

Introduction

Some medical treatments require regular delivery of medications or other fluids into the blood. One method to gain access to the blood is called a vascular access port (also called a port). Ports can be used to remove blood samples for testing and for injection of contrast during certain types of tests known as CT or CECT scans.

This guide provides information about one model of port, SmartPort⁺. Please talk with a healthcare provider if you have questions after reading this information.

What are the SmartPort⁺ implantable ports made of?

The port consists of three parts:

The port body - a small, hollow chamber made of titanium or plastic with a silicone disk (or septum) at the surface for placing a needle. The septum must be accessed with a special needle (non-coring needle) that will allow closing of the septum after the needle is removed. The port has a special feature (Vortex^{*} technology) that allows fluid to reach all surfaces inside it. This helps reduce build up inside the port.

The catheter - a long, soft flexible tube made of polyurethane. One end of the catheter is securely connected to the port, and the other end is placed into a large vein in your chest. Endexo^{*} is a special material that is added to polyurethane to reduce blood component build-up on the surface of the catheter.

Blood clotting reduction is supported by laboratory testing. Results in the lab may be different from the results when used in humans.

The catheter lock - a plastic or metal piece that is used to connect the catheter to the port body.

Materials the ports are made of:

For a complete listing of materials these ports are made of, please visit our website ifu.angiodynamics.com. In the field named "UPN" on the webpage, enter the product part number (a number that begins with "H") that can be found on the Implant Card provided to you when the port was implanted. Press the Search button and locate the Material Information document for your implanted port.

Your healthcare team should know of any drug or material sensitivities you may have prior to port placement and usage.

How is the port placed?


Placement of the port is a surgical procedure. Medications to make the port placement more comfortable, such as sedatives, numbing medications, or pain medications may be given. The medications given during surgery should be discussed with the healthcare team. The skin on the chest (or other location as recommended by the healthcare provider) is cleaned to remove germs. A small incision is made to form a space or “pocket” for the port and a “tunnel” for the catheter.

The port is placed under the skin and into the pocket. One end of the catheter is placed through the tunnel and attached to the port body. The other end of the catheter is placed into a large vein in the chest.

The tip of the catheter is placed near the heart. The tip of the catheter will be checked to make sure the catheter is in the best place before it is used for treatment. The incision is closed with sutures or other methods as determined by the healthcare provider. A dressing may be placed on the incision.

Chest Port Insertion Site

This part of the catheter is in a tunnel under your skin.

This is your chest port site which will be covered with a dressing for a few days until your incision heals.

What to expect during recovery from port placement

Healing may take a few weeks but may be different for each patient. During this time, the area of the surgery may be tender, or sore. This should decrease over time as the body heals.

It is important to follow healthcare providers' instructions on the care of the port insertion site, dressing/bandage, follow up care, and physical activity. The healthcare provider will also decide when the port may be used for treatment.

How is the port accessed and used?

Specially trained healthcare providers can access ports using a special needle (called a “non-coring needle”). These needles deliver medications or fluids and/or can take blood samples. Special care is taken to prepare the skin (sometimes called aseptic or sterile technique) before the needle is placed. This is to reduce bacteria, or germs. The needle is placed through the skin and into the septum of the port. A pricking sensation may be felt when the port is initially accessed. This may become milder over time.

When the port is not in use, it should be flushed from time to time to keep it working normally. Instructions on the routine flushing, care and maintenance schedule should be provided by the healthcare team.

How to identify the SmartPort⁺ device

Ports can be identified in a few ways. Each port kit contains an Implant Card and reminder band that can be worn. This device information should be presented to all healthcare providers, so they are aware of the port. This information is also needed for safe and proper use of the port. It is important to have this information available in case medical care is needed. If the port information is lost, the healthcare provider can provide replacement information.

What is a “CT” or “CECT” Scan, and how is the port used during these types of tests

Contrast Enhanced Computed Tomography (CECT) scans tests provide information about a medical diagnosis and help monitor a medical condition. These tests are sometimes called “CAT” scans or “power injection studies”.

These tests may require injection of contrast material. Before performing this test using an implanted port, it is important that the radiology department confirms which port is in place. SmartPort⁺ may be used for this purpose as long as the manufacturers’ instructions are followed by the healthcare team.

Once the healthcare professional has confirmed that the port can be used for injection of contrast, the port is accessed with a special non coring needle. The needle must also have the capability to inject contrast fluid into blood at the rate needed to perform the test.

How long can a port remain in place?

A port usually remains in place until it is no longer needed for medical care. It can then be removed during a brief surgical procedure.

If the SmartPort⁺ is left in place, AngioDynamics recommends X-rays from time to time. These are done to see if the port is still in the correct location in the body. It can also help find any problems such as pinching of the catheter which may prevent the port from working properly.

Commonly asked questions

Q: How do I take care of my port?

A: Follow the instructions given by your healthcare provider. Be sure to ask if you have questions.

Q: Are AngioDynamics ports latex free?

A: There is no detectable level of natural latex rubber in the SmartPort⁺, associated packaging and contents. These products are not manufactured with natural rubber latex.

Q: Will my port affect my daily activities?

A: Ask your healthcare provider about specific activities that are approved for you, and the safest time to resume them.

Q: Will I need to wear a bandage over my port?

A: It is important to follow healthcare providers' instructions on the care of the port insertion site, dressing/bandage, follow up care, and physical activity. The healthcare provider will also determine when the port may be used for treatment.

Q: Will my port activate security alarms?

A: Security systems will most likely not detect the small amount of metal in the device. If needed, simply show your patient identification card. You may declare your implant prior to entering security screening in order to be screened in the most appropriate way. If you do not have a Patient Implant card, please call 1-800-772-6446.

Q: Is SmartPort⁺ port safe with CECTs?

A: The materials used in the SmartPort⁺ devices are safe with CECT procedures up to 300 psi when used according to the manufacturer's instructions.

Q: Is SmartPort⁺ safe for MRIs studies?

A: MRI is safe up to 7T. Please ask your healthcare provider if you have any questions.

Q: What if my healthcare provider has not seen a patient with a SmartPort⁺ device before?

A: Always show healthcare providers your Patient Implant Card as it contains a summary of important information about your port. If they have more questions, they can call the AngioDynamics Clinical Information number at 1-800-772-6446 option 5.

Q: Are non-coring infusion sets latex free?

A: The healthcare team should refer to the packaging and label of the non-coring needle to determine if it contains latex rubber.

Q: How often should the non-coring needle be changed if the patient is receiving ongoing therapy?

A: The non-coring needle may be changed every 7 days, or as directed by the healthcare provider.

Important information to report to your healthcare provider

It is important to talk with a healthcare provider about any illnesses you may have before getting a port or after the procedure. A list of potential problems is provided with the device. Talk to your healthcare team or call AngioDynamics to know more about this information.

Risks and benefits of the getting the port should be discussed with your healthcare team.

It is important to notify your healthcare provider(s) right away about the following, as well as any new symptoms that you are concerned about:

- Aching, swelling or pain in the shoulder, neck, or arm.
- Pain, redness, swelling, or soreness over or around the port site after the incision heals.
- Any unusual drainage from the incision site.
- Burning sensation upon infusion of medication or fluids.
- Chills, shortness of breath, dizziness, or fever.

These may be signs of a complication related to the port or it's placement.

Patient checklist

- ☐ Keep your Patient Implant Card with you.
- ☐ Present the Patient Implant Card to any healthcare providers involved in your care.
- ☐ If the healthcare provider needs additional information, they may call AngioDynamics at 1-800-772-6446.
- ☐ Discuss any concerns or questions about your port, or this patient guide, with your healthcare provider.
- ☐ Know your treatment schedule - when, where, and how often you will go for treatment.
- ☐ If traveling, discuss port care and treatment for your trip.
- ☐ Follow healthcare providers' instructions on the care of the port insertion site, dressing/bandage, follow up care, and physical activity. A healthcare provider will decide when the port may be used for treatment.
- ☐ Additional patient information and educational tools may be accessed by visiting www.angiodynamics.com or by calling AngioDynamics at 1-800-772-6446.

Important Notes and Instructions from the Healthcare Provider:

Surgical Site and Dressing Care Instructions:

[illegible]

Post Op Physical Activity Instructions:

[illegible]Post Op Follow Up Appointment Day, Time, and Location:

Date Port Cleared for Use for Treatment or Tests:

Routine Care and Maintenance Schedule:

Date, Location and Time of First Treatment:

Treatment Schedule:

Date, Location and Time of Blood Tests:

Date, Location and Time of Other Tests:

Healthcare Provider Contact Information:

Emergency Contact Information:

Commonly used terms

Implantable Port:

A fluid or medication delivery device with a hollow chamber, called the port. The port is connected to a soft, hollow tube called the catheter.

Power Injectable Port:

A type of implantable port that provides access for both IV therapy and for CECT tests.

Contrast-Enhanced

Computed Tomography (CECT):

A type of x-ray test that requires a pump to deliver testing fluid fast and at a high pressure. These tests produce superior images of your body to help your medical team better manage your treatment.

Septum:

The septum is the center of the port. A non-coring needle is inserted through the septum to deliver your fluid or medications.

Catheter:

Hollow tube connected to the port body that is placed into a blood vessel.

Infusion:

The delivery of fluid and medication into a blood vessel.

Magnetic Resonance Imaging (MRI):

A type of scan that uses strong magnetics and radio waves to create detailed pictures of the inside of the body.

Symbols for US: In compliance with the requirements of 21CFR Part 801.15, a glossary of symbols which appear without accompanying text within the product labeling is provided below

Symbol	Ref	Title of Symbol	Meaning of Symbol
	5.7.7	Medical device	Indicates the item is a medical device. ^a
	5.7.10	Unique Device Identifier	Indicates a carrier that contains unique device identifier information. ^a
	3.1.11	Magnetic Resonance Conditional	An item with demonstrated safety in the MR environment within defined conditions including conditions for the static magnetic field, the time-varying gradient magnetic fields and the radiofrequency fields. ^b
	NA	Contrast Enhanced Computed Tomography	Contrast Enhanced Computed Tomography

a. EN ISO 15223-1 - Medical Devices – Symbols to be used with medical device labels, labeling and information to be supplied.

b. ASTM F2503-20 - Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment

Notice for patients in Australia:

Please report product concerns/complaints to the Therapeutic Goods Administration (TGA) at <https://www.tga.gov.au/reporting-adverse-events>.

Notification of adverse event or complaint to AngioDynamics or the regulatory body in your country is not a replacement for medical opinion

SMARTPORT⁺

INDICATIONS FOR USE:

The ports are indicated for patients who require long-term access to the central venous system for blood specimen withdrawal and administration of fluids including but not limited to hydration fluids, chemotherapy, analgesics, nutritional therapy and blood products, as well as the administration and adequate removal of nuclear medicine.

When used with power injectable needles, the ports are indicated for power injection of contrast media. For power injection of contrast media, the maximum recommended infusion rate is 5 mL/s with 19G or 20G non-coring power injectable needles or 2 mL/s with a 22G non-coring power injectable needle.

Needle Size (G), non-coring power injectable	Catheter Size (F)	Maximum Recommended Flow Rate Setting (mL/s)	Maximum Recommended Pressure Setting (psi)
19/20	5, 6, and 8	5	300
22	5, 6, and 8	2	300

Considering the indication for long-term access to the vascular system, the ports may, at the discretion of the clinician, be employed in the clinical management, evaluation, and monitoring of conditions such as, but not limited to, the following:

- Dehydration
- Malnutrition
- Malignant conditions
- Fluid and electrolyte imbalances
- Pain syndromes
- Anemia
- Hematologic abnormalities

CAUTION: Federal Law (USA) restricts this device for sale by or on the order of a physician

NOTE TO HEALTHCARE PROFESSIONAL:

Refer to manufacturer's Instructions for Use for full prescribing and safety information

- Always use a sterile technique for skin preparation
- Always use a non-coring needle to access the port
- After use, flush per institutional protocol
- DO NOT exceed a 300 psi pressure limit setting or the maximum flow rate setting indicated on the power injection machine if power injecting through the port.

* AngioDynamics, the AngioDynamics logo, SmartPort, the SmartPort logo, and Vortex are trademarks and/or registered trademarks of AngioDynamics, Inc., an affiliate or a subsidiary. Other trademarks are property of their respective owners.

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All rights reserved.
© 2025 AngioDynamics, Inc., or its affiliates.