angiodynamics

# AngioVac Cannula

|                           |    |
|---------------------------|----|
| Directions For Use .....  | 3  |
| Instrucciones de uso..... | 10 |
| Mode d'emploi.....        | 17 |
| Gebrauchsanweisung .....  | 24 |
| Istruzioni per l'uso..... | 31 |
| Gebruiksaanwijzing .....  | 38 |
| Brugsanvisning.....       | 45 |

16600506-01

A

2020-01

| <b>TABLE OF CONTENTS</b>   | <b>ENGLISH</b> |
|--|----------------|
| <b>WARNING</b> .....   | <b>3</b>       |
| <b>DEVICE DESCRIPTION</b> .....                                    | <b>3</b>       |
| <b>INDICATIONS FOR USE</b> .....                                   | <b>3</b>       |
| <b>INTENDED USE</b> .....  | <b>3</b>       |
| <b>CONTRAINDICATIONS</b> .....                                     | <b>3</b>       |
| <b>WARNINGS</b> .....  | <b>4</b>       |
| <b>PRECAUTIONS</b> .....   | <b>4</b>       |
| <b>ADVERSE EVENTS</b> .....  | <b>5</b>       |
| <b>HOW SUPPLIED</b> .....  | <b>5</b>       |
| Graph: Pressure Drop vs Flow Rate.....                             | <b>5</b>       |
| <b>OPERATIONAL INSTRUCTIONS</b> .....                              | <b>5</b>       |
| <b>TROUBLESHOOTING SUSTAINED LOW FLOW/NO FLOW CONDITIONS</b> ..... | <b>7</b>       |
| Specification Table .....  | <b>7</b>       |
| <b>WARRANTY</b> .....  | <b>8</b>       |

# AngioVac Cannula

## Rx ONLY

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

---

### WARNING

**Contents supplied STERILE using an ethylene oxide (EO) process.** Do not use if sterile barrier is damaged. If damage is found, call your sales representative. Inspect prior to use to verify that no damage has occurred during shipping.

**For single patient use only. Do not reuse, reprocess or resterilize.** Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

---

### DEVICE DESCRIPTION

The *AngioVac*\* Cannula has a sheath deployed self-expandable, nitinol reinforced, funnel-shaped distal tip to facilitate flow, and to facilitate en bloc removal of undesirable intravascular material, including, but not limited to, thrombus, embolus, clot or debris, during the performance of venovenous extracorporeal circulation. Target vessels for thrombus/ embolus extraction include, but are not limited to, the iliofemoral vein, Inferior Vena Cava (IVC), Superior Vena Cava (SVC) and Right Atrium (RA). The nitinol basket and the radiopaque markers located in the distal tip of the funnel and sheath aid in device visualization. Standard surgical or percutaneous insertion techniques can be employed. If clinically indicated, the cannula inner diameter can accommodate devices (e.g., angiographic catheters) with a maximum outer diameter of 17F (5.7mm/0.223in.). The cannula is available in angled configurations of ~20° and ~180°.

### INDICATIONS FOR USE

The *AngioVac* Cannula is indicated for use as a venous drainage cannula and for removal of fresh, soft thrombi or emboli during extracorporeal bypass for up to 6 hours.

### INTENDED USE

The *AngioVac* Cannula is intended to be used with the *AngioVac* Circuit and commonly available extracorporeal equipment (e.g., centrifugal pump, bubble trap (filter) and reinfusion cannula) to facilitate en bloc removal of undesirable intravascular material during venovenous extracorporeal circulation.

### CONTRAINDICATIONS

The following contraindications are applicable:

- Do not use if the patient has severe arterial or venous vascular disease.
- The device is contraindicated in the removal of chronic firmly adherent intravascular material (e.g., atherosclerotic plaque, chronic pulmonary embolism).
- The device is contraindicated for use in the right heart or pulmonary arteries during active cardiopulmonary resuscitation.

---

## WARNINGS

- Inspect product prior to, during, and after use, to ensure that no damage to the product has occurred. Failure to inspect could result in injury to patient or user.
- Verify that the lumen is patent and that the cannula has not been damaged or kinked prior to use.
- Instructions for use and manuals for the AngioVac Cannula and all related accessories and extracorporeal circulatory devices should be read prior to use and devices used as indicated.
- Selection of the patient as a candidate for use with this device and for such procedures as it is intended is the physicians' sole responsibility. The outcome is dependent on many variables including, patient pathology, surgical procedure, and perfusion procedure/technique. The benefits of use of this device must be weighed against the risks including risks of systemic anticoagulation and must be assessed by the prescribing physician.
- As with all medical devices, this device and ancillary equipment are to be used by trained physicians only. Specifically, this device is to be used only by medical personnel experienced with initiating and monitoring extracorporeal bypass and by physicians trained and experienced using surgical and/or percutaneous (Seldinger) vascular access techniques.
- User assumes responsibility for any variations in the extracorporeal circulation/procedures that could compromise the intended use of this device.
- Do not use in conjunction with a power injector.
- Do not alter the AngioVac Cannula in any way.

---

## PRECAUTIONS

- The AngioVac Cannula can only be used in conjunction with the AngioVac Circuit.
- A strict anticoagulation protocol should be followed and anticoagulation should be carefully monitored during the procedure.
- Prior to use, confirm that the cannula size is appropriate for the vessel to be accessed and for all instruments, catheters, sheaths and cannula used during the procedure.
- Cannula tip placement and positioning should be guided and confirmed using standard guide wire and fluoroscopic techniques.
- Do not insert, attach or manipulate cannula in a manner that may result in extravascular placement, kinking or compression or in a way that may alter or restrict flow.
- Use caution when connecting the cannula to the AngioVac Circuit so as not to damage the cannula or circuit.
- **DO NOT** clamp the cannula. Clamping of the cannula may result in permanent cannula wall distortion and/or lumen collapse.
- Ensure that, once inserted and in an appropriate intravascular position, the cannula is completely free of air by aspirating using the luer lock port on the cannula sheath prior to connecting to the extracorporeal circuit.
- Caution should be used in positioning the cannula as undue pressure exerted by the tip can cause perforation of, or damage to vessels and intravascular structures. If difficulty or resistance is encountered during placement or withdrawal, the cause should be determined and corrected before proceeding. Failure to do so may result in damage to the vessel.
- When positioning the cannula avoid impingement of the tip against vessel walls or vessel side branches as this may obstruct flow and/or damage the vessel.
- To minimize the potential for vascular trauma ensure that the cannula is placed in an appropriately sized vessel.
- Do not attempt to insert a catheter or other device having a diameter larger than that which is compatible with the cannula (i.e. 17F (5.7 mm/0.223 in). Device/cannula damage or breakage may occur.
- The use of cold solutions administered externally may increase the stiffness of the cannula and alter/increase the pressure exerted on the tip during manipulation.
- **DO NOT** use alcohol or alcohol-based fluids for lubrication as these solutions may damage the cannula components.
- During use, carefully monitor cannula for both inflow and outflow obstruction/occlusion.

## ADVERSE EVENTS

This device, as do all extracorporeal blood vessel devices, has possible side effects, which include but are not limited to infections, blood loss, thrombus formation, embolic events, vessel, ventricular or valvular damage and complications of percutaneous or surgical insertion techniques. These may occur if the instructions for use are not followed.


Possible complications include those normally associated with large bore surgical and/or percutaneous vessel catheterization/cannulation, anticoagulation, extracorporeal circulation and application of intravascular introducer systems which include but are not limited to:

- Air Embolism
- Arrhythmias
- Blood Loss/Blood Trauma
- Death
- Damage to Vessel
- Device Fracture
- Distal Embolization of Thrombus
- Hematoma
- Hemoptysis
- Hemorrhage
- Local or Systemic Infections
- Pleural Effusion
- Pulmonary Embolism
- Pulmonary Infarction
- Vascular Thrombosis
- Venous Valve Injury
- Ventricular Perforation
- Vessel Spasm

## HOW SUPPLIED

Contents supplied STERILE using an Ethylene Oxide (EO) Process. Store in a cool, dry, place. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

### Graph: Pressure Drop vs Flow Rate

## OPERATIONAL INSTRUCTIONS

1. The patient should be prepared and draped in the usual and typical sterile manner for percutaneous/surgical vascular procedures.
2. Sterile vascular cannulation should be performed utilizing typical percutaneous or open surgical techniques.
3. Carefully open the sterile pouch containing the cannula and obturator and remove from their protective packaging without damaging the device.
4. Assemble the AngioVac Cannula.
  - a) Collapse the AngioVac funnel tip by pulling the Cannula into the sheath
  - b) Tighten the tuohy on the sheath hub to secure the cannula in the collapsed position
  - c) Connect cannula to the primed "Y" AngioVac circuit hub.
  - d) Hydrate the obturator with saline solution and flush ID with syringe
  - e) Remove the wire tuohy insert tool and fully insert the obturator into the red tuohy and tighten tuohy until finger tight.

- f) Flush the AngioVac Sheath
  - Connect a syringe filled with saline to the stopcock on the sheath.
  - Turn the stopcock to allow fluid flow from the syringe to the sheath.
  - Flush the sheath to dispel all air between the cannula and the sheath.
  - Turn the stopcock to close the fluid path to the sheath.
- g) Open perfusionist clamp at drainage side of hub to flush cannula. Close clamp.
5. If desired for use, assemble the introducer sheath/introducer dilator.
  - a) Lubricate the introducer dilator with saline solution or appropriate water soluble lubricant.
  - b) Insert the introducer dilator into the introducer sheath until the hub is well seated and secured in that position.
6. Identify the vessel to be cannulated and ensure that it is of adequate size to allow for introduction of the introducer sheath and/or AngioVac Cannula. In a percutaneous fashion or into the exposed vessel, introduce an 18g needle into the vessel and aspirate blood. If free flow of blood is not present upon aspiration, reposition the needle and aspirate. Repeat as needed until adequate blood flow is established.
7. Insert guidewire through the needle into the vessel.

---

**CAUTION:** The guidewire should advance without resistance and should be done under fluoroscopic guidance.

---

**CAUTION:** Do not withdraw the guidewire back into the needle in the vessel as shearing of the wire may occur. Always remove the needle first.

---

8. Advance the guidewire into the vessel and position it at or just beyond the desired cannula tip position.
9. Hold the guidewire in place and remove the introducer needle by sliding it away from and out of the patient over the guidewire.
10. Enlarge the insertion site of the wire into the skin or vessel with an #11 blade scalpel to approximate the size of the sheath and/or cannula.
11. Hold the guidewire firmly in position and pass a smaller stepped dilator over the guidewire and into the vessel to enlarge the tract to the largest part of the dilator
12. Remove the vessel dilator by sliding it away from and out of the patient over the guidewire. Be sure to apply digital pressure to the entrance site for hemostasis and hold the guidewire firmly in position.
13. Repeat steps 11 and 12 with larger dilators as appropriate to allow placement of the introducer sheath/AngioVac Cannula.
14. If an introducer sheath is used pass the introducer sheath/dilator assembly over the guidewire, advancing until the guidewire can be grasped beyond the introducer hub.
15. Holding the guidewire firmly in position, advance the introducer sheath/dilator assembly over the guidewire into the vessel to the desired depth.

---

**CAUTION:** Advancement of the introducer sheath/dilator assembly with a rotating motion make placement easier. Do not advance if resistance is met as doing so may result in product and/or vessel damage.

---

16. While holding the guidewire and introducer sheath in place, remove the dilator from the sheath sliding it over the guidewire away from and out of the patient.
17. If applicable, secure the introducer sheath in place by suturing it to the skin at the point of insertion.
18. Place the AngioVac Cannula/obturator over the guidewire and advance until the guidewire can be grasped as it exits the distal end of the obturator.
19. While holding the guidewire firmly in position, advance the AngioVac Cannula/obturator over the guidewire and through the introducer sheath into the vessel.
20. If the AngioVac Cannula is used without an introducer sheath, once the vessel and/or insertion tract has been dilated to the appropriate size, place the AngioVac Cannula/obturator assembly over the guidewire and advance into the vessel while holding the guidewire firmly in the appropriate position.

---

**CAUTION:** Advancement of the AngioVac Cannula/obturator assembly with a rotating motion makes placement easier. Do not advance if resistance is met as doing so may result in product and/or vessel damage.

---

21. Using fluoroscopy, advance the AngioVac Cannula/obturator into the desired position over the guidewire.
22. Once the optimal cannula position has been achieved, remove the guidewire by sliding it out of and away from the patient through the obturator.
23. Once the guidewire has been removed, remove the obturator from the AngioVac Cannula completely by loosening red tuohy and pulling it out and away from the patient.
24. Insert wire tuohy adapter and tighten red tuohy until finger tight.

---

**CAUTION:** Perfusion tubing clamps must be used. Other clamps may be too weak to ensure complete occlusion of the Y-Connection. **DO NOT** clamp the wire reinforced cannula.

---

25. Connect a syringe to the side port on the cannula. Aspirate air from the space between the sheath and cannula. Then re-flush. (If air found, remove via other port on stopcock.) Remove syringe and close valve on stop cock.
26. Loosen sheath tuohy to allow cannula movement. Advance the AngioVac Cannula out of the sheath expanding the funnel tip and exposing the desired tip angle.
27. Tighten the tuohy on the sheath onto the cannula locking the desired tip angle in place (adjust placement as needed). Ensure all appropriate clamps are again in the appropriate open or closed position and that there is no inflow or outflow obstruction.
28. Once adequate anticoagulation has been achieved, initiate extracorporeal circulation as desired. Engage undesirable intravascular material (UIM) to distal end of cannula.
29. Pull back cannula tip into sheath prior to removal of cannula assembly. Remove AngioVac Cannula Assembly from patient.

#### **TROUBLESHOOTING SUSTAINED LOW FLOW/NO FLOW CONDITIONS**

A sudden decrease in flow may be a signal that a large mass has become seated in the funnel. Once the mass is compressed and extracted through the cannula flow will be re-established. If flow is not quickly re-established, the following maneuvers should be performed:

1. Reposition the cannula away from the clot to an anatomic position where good flow was previously achieved to see if flow can be re-established.
2. If the low flow/no flow condition is sustained, **WITH THE PUMP OPERATIONAL**, pull back cannula tip into sheath and remove the device from the venous entry site to check for luminal obstruction. The cannula may be flushed by either disconnecting the cannula from the circuit and flushing with saline using a bulb syringe, or keeping the circuit connected clamp the priming line and apply pressure on the IV prime bag to flush the cannula.

---

**CAUTION:** Assure that all air is removed from circuit/cannula prior to re-establishing venous access.

---

#### **Specification Table**

|                                    |                    |
|------------------------------------|--------------------|
| Cannula Working Length             | 77 cm              |
| Recommended Introducer Sheath Size | 26 Fr              |
| Recommended Guidewire              | 0.038 in (0.97 mm) |
| Approx. Priming Volume             | 30 mL              |
| Funnel Tip Open Diameter           | 14 mm              |

## WARRANTY

AngioDynamics warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond AngioDynamics control directly affect the instrument and the results obtained from its use. AngioDynamics obligation under this warranty is limited to the repair or replacement of this instrument and AngioDynamics shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. AngioDynamics neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **AngioDynamics assumes no liability with respect to instruments reused, reprocessed, resterilized, modified or altered in any way and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

\* AngioDynamics, the AngioDynamics logo, AngioVac, and the AngioVac logo are trademarks and/or registered trademarks of AngioDynamics, Inc., an affiliate or a subsidiary.

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© 2020 AngioDynamics, Inc. or its affiliates. All rights reserved.