# Solero

## Cart Owner's Manual

Directions For Use.....	2
Instrucciones de uso.....	4
Mode d'emploi .....	6
Gebrauchsanweisung.....	8
Istruzioni per l'uso.....	10
Gebruiksaanwijzing .....	12
Brugsanvisning .....	14
Bruksanvisning .....	16
Οδηγίες χρήσης.....	18
Instruções de utilização .....	20
Használati utasítás .....	22
Pokyny k použití .....	24
Wskazania do stosowania .....	26
Bruksanvisning .....	28
Kullanım Yönergeleri .....	30
Naudojimo instrukcijos.....	32
Instrucțiuni de utilizare.....	34
Показания к применению .....	36
Käyttöohjeet.....	38
Български. ....	40

# Solero

## Cart Owner's Manual

ENGLISH

**⚠️ ONLY Caution:** Federal Law (USA) restricts this device to sale by or on the order of a Physician.

### DEVICE DESCRIPTION:

The *Solero*\* MTA Cart is an accessory device to the Solero Microwave Tissue Ablation (MTA) System. It is intended for transport and operation of the Solero Generator and is used in control environments such as those found within surgical suites and emergency rooms.

### INDICATIONS FOR USE:

In the United States the Solero MTA System has the following Indications for Use:

- The Solero Microwave Tissue Ablation (MTA) System and accessories are indicated for the ablation of soft tissue during open procedures. The Solero MTA System is not intended for cardiac use.

The Solero MTA Cart is CE marked as an accessory device to the Solero Microwave Tissue Ablation (MTA) System.

### CAUTIONS AND WARNINGS

- Read the Operator's Manual and Directions for Use for the Solero MTA System and Accessories prior to use. Safe tissue ablation is dependent upon operation of the equipment in accordance with the operating instructions.
- Discontinue use and contact AngioDynamics if the cart is damaged.
- Do not put more than 5 kg/11 lbs. weight in tray.
- Do not hang more than 3 Liters of saline on the IV Pole.
- Do not push or transport cart by gripping the IV pole or the Solero MTA Generator.
- Do not use for storage and transport of other medical devices.
- Do not sterilize the Solero Cart.
- Do not use if package is opened or damaged.
- Do not use if labeling is incomplete or illegible.
- Ensure that the Solero MTA Generator is properly seated on the cart.
- Lock wheels when not in transport. Unlock wheels when transporting.
- The Solero MTA System and accessories are MR Unsafe.

Any serious incident which has occurred with the use of this device should be reported to AngioDynamics at [complaints@angiodynamics.com](mailto:complaints@angiodynamics.com) and to the National Competent Authority. Refer to the following web address for contact information of the EU Competent Authority. [https://ec.europa.eu/health/sites/health/files/md\\_sector/docs/md\\_vigilance\\_contact\\_points.pdf](https://ec.europa.eu/health/sites/health/files/md_sector/docs/md_vigilance_contact_points.pdf) The instructions for use are available electronically at [ifu.angiodynamics.com](http://ifu.angiodynamics.com)

### CART SPECIFICATIONS:

Physical Dimensions (Nominal):	
Width:	493 mm/19.4 in
Depth:	557 mm/21.9 in
Height to Shelf:	889 mm/35.0 in
Height to IV Pole:	1706 mm/67.2 in
Weight:	19 kg /42 lbs
Max weight of Cart, Solero Generator, Tray Contents, and Saline Bags is approximately 45 kg/99 lbs.	

### SYSTEM SET UP AND PREPARATION

#### Receipt:

- Remove the Solero MTA Cart and Owner's Manual from the shipping container.
- Visually inspect the contents for damage during shipment.
- Thoroughly read the Owner's Manual. For Assembly instructions refer to page 29.

#### Transportation:

- To transport the system within the hospital, first ensure that the power cable to the Solero MTA Generator is disconnected from the supply mains electrical outlet and the wheels are unlocked.
- Ensure that all cables are properly stored on the cart.
- Move the cart by the handle. Avoid gripping the IV pole or the Solero MTA Generator.

#### Cleaning:

- Clean the cart by wiping down all accessible surfaces with a clean cloth that is slightly dampened with 70% isopropyl alcohol cleaning solution.

#### User Replacement Part Kits:

- The Solero MTA Cart is user serviceable. Replacement part kits can be purchased through AngioDynamics, see back page for contact information.

Description	Part Number	Assembly Diagram
IV Pole	12740140	Page 30
Power Cord Hook	12740141	Page 31
Storage Tray	12740142	Page 32
Caster	12740143	Page 33

### WARRANTY

#### Limited Warranty For Angiodynamics Products

This limited warranty applies to AngioDynamics, Inc. ("AngioDynamics") products including its Solero Generators, carts, and accessories ("Products"). "Purchaser" as used herein refers to any individual or entity that purchases any Product from AngioDynamics or its authorized representative.

- AngioDynamics agrees to repair or replace (in AngioDynamics' sole discretion) any Product that is proven to have a material defect in material or workmanship and of which AngioDynamics is notified in writing prior to the limited warranty expiration date, if applicable, or prior to the date indicated on the warranty card accompanying the Product. If the Product does not have an expiration date or if no warranty card is available, AngioDynamics' obligation to repair or replace the Product shall not exceed 12 months from the receipt of such Product by the individual or entity originally purchasing such Product directly from AngioDynamics or directly from AngioDynamics' authorized representative. Purchaser's sole and exclusive remedy against AngioDynamics, and AngioDynamics' sole and exclusive liability, shall be the repair or replacement of the Product in accordance with this limited warranty.
- Except as set forth below, the Products contain no user-serviceable parts. If servicing is required, the Products must be returned to AngioDynamics and may only be returned with the prior written approval of AngioDynamics. Any such approval shall reference a Service Order number issued by AngioDynamics Customer Service. Shipping and transportation costs, if any, incurred in connection with the return of defective Product to AngioDynamics shall be the responsibility of the Purchaser. If the Product is determined, in AngioDynamics' sole discretion, to be defective, shipping and transportation costs will be reimbursed by AngioDynamics to the Purchaser.
  - The following user-serviceable parts are available from AngioDynamics for the Solero cart: (a) storage tray; (b) caster; (c) power cable hook; and (d) IV pole. Purchaser can order these parts from AngioDynamics and replace the parts itself without voiding this limited warranty.
- Except for the express limited warranties set forth in Section 1 above, AngioDynamics grants no warranties for Products, express or implied, either in fact or by operation of law, by statute or otherwise, and AngioDynamics specifically disclaims any implied warranty of quality, warranty of merchantability, warranty of fitness for a particular purpose or warranty of non-infringement.
- AngioDynamics' liability under these warranties shall be limited to a refund of the Purchaser's purchase price or repair or replacement of the Product. In no event shall AngioDynamics be liable for the cost of procurement of substitute goods by the Purchaser or for any special, consequential or incidental damages for breach of warranty.
- AngioDynamics shall not be liable, expressly or implied, for:
  - Repairs or modifications to the Products performed other than by AngioDynamics or an AngioDynamics authorized repair facility, except as expressly set forth in this limited warranty; or

5.2. Use of the Products in any manner or medical procedure, other than that for which it is designed, either of which shall void this warranty.

\* AngioDynamics, the AngioDynamics Logo, Solero and the Solero logo are trademarks and / or registered trademarks of AngioDynamics, Inc., an affiliate or a subsidiary.

In compliance with the requirements of 21 CFR Part 801.15, below is a glossary of symbols which appear without accompanying text.

Symbol	Ref	Title of Symbol	Meaning of Symbol
	5.1.1	Manufacturer	Indicates the medical device manufacturer. <sup>a</sup>
	5.1.2	Authorized representative in the European Community	Indicates the authorized representative in the European Community. <sup>a</sup>
	5.1.3	Date of manufacture	Indicates the date when the medical device was manufactured. <sup>a</sup>
	5.1.7	Serial number	Indicates the manufacturer's serial number so that the medical device can be identified. <sup>a</sup>
	5.1.8	Importer	Indicates the entity importing the medical device into the locale. <sup>a</sup>
	5.3.1	Fragile, handle with care	Indicates a medical device that can be broken or damaged if not handled carefully. <sup>a</sup>
	5.3.4	Keep dry	Indicates a medical device that needs to be protected from moisture. <sup>a</sup>
	5.7.7	Medical device	Indicates the items is a medical device. <sup>a</sup>
	5.7.10	Unique Device Identifier	Indicates a carrier that contains unique device identifier information. <sup>a</sup>
<b>UPN</b>	NA	Universal Product Number	A Universal Product Number (UPN) code represents the manufacturer's number for an item.
<b>CE</b>	NA	CE Mark	Manufacturer's declaration of conformity to the Medical Device Regulation EU 2017/745. <sup>b</sup>
	NA	Magnetic Resonance (MR) unsafe	Keep away from magnetic resonance imaging (MRI) equipment. <sup>c</sup>

Symbol	Ref	Title of Symbol	Meaning of Symbol
	M002	Follow instructions for use ifu.angiodynamics.com	Refer to Instruction Manual or booklet. <sup>d</sup>
	P017	No pushing	To prohibit pushing against an object. <sup>d</sup>
	1135	Recyclable Package	Recyclable Package. <sup>e</sup>
	1321A	Mass, weight	To indicate mass. <sup>f</sup>
	0623	This way up	This is the correct upright position of the distribution packages for transport and/or storage. <sup>f, g</sup>
	2402	Do not stack	Stacking of the distribution packages is not allowed and no load shall be placed on the distribution packages. <sup>f, g</sup>

- a. EN ISO 15223-1 - Medical Devices – Symbols to be used with medical device labels, labeling and information to be supplied.
- b. EU 2017/745 Medical Device Regulations published 5 May 2017
- c. ASTM F2503-13 - Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.
- d. BS EN ISO 7010 Graphical Symbols
- e. EN ISO 14021 Environmental labels and declarations. Self-declared environmental claims (Type II environmental labeling)
- f. ISO 7000 - Graphical symbols for use on equipment – Registered symbols.
- g. ISO 780 Distribution packaging - Graphical symbols for handling and storage of packages

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AngioDynamics UK Ltd.