Uni-Fuse INFUSION SYSTEM

CE 2797

INSTRUCTIONS FOR USE

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2022-04

INTENDED USE

The Uni-Fuse[®] Infusion System is intended for the administration of thrombolytic agents into the peripheral vasculature.

INDICATIONS

The Uni-Fuse[®] Infusion System is intended to be used to treat peripheral venous thrombosis and peripheral arterial thrombosis by catheter directed thrombolysis.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged or unintentionally opened before use. If damage is found, call your sales representative. Inspect prior to use to verify that no damage has occurred during shipping.

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Uni-Fuse Infusion System is to be treated as contaminated biomedical waste subsequent to use. The used or unused devices, should be disposed of in accordance with hospital, administrative and/or local government policy for such items.

Untampered device packaging should be recycled if applicable, or disposed of as common waste in accordance with hospital, administrative and/or local government policy for such items.

PRODUCT DESCRIPTION

The Uni-Fuse Infusion System consists of the following components:

- PRO[®] (Pressure Responsive Outlet) Infusion Catheter with longitudinal slits at the distal end. Radiopaque markers on the catheter shaft indicate the active infusion pattern.
- Occluding Ball Guidewire which occludes the distal end of the PRO Infusion Catheter.

The infusion patterns are determined by the catheter shaft slit lengths. The 2cm, 5cm, 10cm, 20cm and 30cm slit lengths include slits at 90° intervals around the axis of the catheter shaft and rows every 5mm. The 40cm and 50cm slit lengths include two slits opposed at 180° placed every 10mm along the catheter shaft.

Introducer Sheath compatibility:

- 4F catheter shaft max. OD 0.058 inches (1.48mm) is compatible with Introducer Sheath ID > 0.058 inches (1.48mm)
- 5F catheter shaft max. OD 0.069 inches (1.76mm) is compatible with Introducer Sheath ID > 0.069 inches (1.76mm)

CONTRAINDICATIONS

The Uni-Fuse Infusion System is contraindicated for use in the coronary and cerebral vasculature.

The Uni-Fuse Infusion System is not intended for the infusion of blood or blood products. Refer to the product insert of the therapeutic solution for indications, contraindications, side effects, cautions and warnings.

POTENTIAL COMPLICATIONS

Potential Complications include, but are not limited to:

- Hematoma at the entry site
- Vessel perforation
- Vasospasm
- Hemorrhage
- Embolism
- Drug reaction
- Neurological deficits including stroke and death
- Pain and tenderness
- Vessel dissection
- Vascular thrombosis
- Other - case (or patient) to which the therapeutic solution is inadequate
- Infection

WARNINGS

- Single patient use. Do not re-use or resterilize.
- Do not use if package is damaged. Do not use if labeling is incomplete or illegible.
- This product should be used only by physicians who have a thorough understanding of angiography and percutaneous interventional procedures.
- AngioDynamics cannot guarantee the proper function of these components if any component is substituted by those of another manufacturer.
- Failure to use an introducer sheath may result in damage to the catheter or occluding ball wire.
- Reprocessing may compromise the integrity of the device and/or lead to device failures.
- Guidewire contains cobalt. Cobalt is classified as CMR 1B and is present in a concentration above 0.1% weight by weight.

Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient.

CAUTIONS

Rx ONLY Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

Intended Patient Population

The Uni-Fuse Infusion System is intended for patients requiring short-term peripheral venous and arterial access for administration of thrombolytic agents. Peripheral venous and arterial access may be required for therapeutic treatment of peripheral venous and arterial thrombus.

Clinical Benefits

The Intended Clinical Benefits of the Uni-Fuse Infusion System include those therapeutic applications afforded by peripheral venous and arterial access. Specifically, the Infusion catheters intended clinical benefits include catheter directed thrombolysis for patients requiring peripheral venous and arterial access for administration of thrombolytic agents resulting in;

- Maximized delivery of the thrombolytic resulting in better vessel patency compared to anticoagulation alone.¹

Serious Incidences

Any serious incident which has occurred with the use of this device should be reported to AngioDynamics or your in-country representative.

Notices for European Union Only:

- Any serious incident which has occurred with the use of this device should be reported to AngioDynamics at complaints@angiodynamics.com and to the National Competent Authority. Refer to the following web address for contact information for the Competent Authorities. EU: https://ec.europa.eu/health/sites/health/files/md_sector/docs/md_vigilance_contact_points.pdf

The physician should discuss with the patient the risks associated with the device. The instructions for use are available electronically at ifu.angiodynamics.com

HOW FURNISHED


Sterile, disposable. Single patient use. Store in a cool, dry place. Maximum storage temperature 27°C.

INSTRUCTIONS FOR USE

1. Use aseptic technique.
2. Flush the PRO Infusion catheter with sterile, heparinized normal saline.
WARNING: Complications may occur, if all the air has not been removed from the infusion catheter and displaced with saline prior to insertion into the body.
3. Insert the PRO catheter to the desired location in the peripheral vascular system using fluoroscopic guidance.
4. Insert the occluding ball wire into the PRO catheter and attach to the PRO catheter. Tighten by hand.
WARNING: Never advance the guidewire against resistance; this could cause vessel trauma and/or wire damage. Determine the cause of resistance under fluoroscopy and take any necessary remedial action.
5. With syringe or slow infusion pump pulse or slow infuse through the proximal luer lock.

Uni-Fuse INFUSION SYSTEM COMPONENTS

- Occluding Ball Guidewire


WARRANTY

AngioDynamics, Inc. warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other representations and warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond AngioDynamics, Inc.'s control directly affect the instrument and the results obtained from its use. AngioDynamics, Inc.'s obligation under this warranty is limited to the repair or replacement of this instrument and AngioDynamics, Inc. shall not be liable for any incidental or consequential losses, damages or expenses directly or indirectly arising from the use of this instrument. AngioDynamics, Inc. neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. AngioDynamics, Inc. **assumes no liability with respect to instruments reused, reprocessed or resterilized, modified or altered in any way, and makes no warranties or representations, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

¹ Elsharawy M, Elzayat E. Early results of thrombolysis vs anticoagulation in iliofemoral venous thrombosis. A randomised clinical trial. *European Journal of Vascular and Endovascular Surgery*. 2002;24(3):209-14.

* AngioDynamics, the AngioDynamics logo, Uni-Fuse and PRO are trademarks and/or registered trademarks of AngioDynamics, Inc., an affiliate or a subsidiary.

Symbols for US: In compliance with the requirements of 21CFR Part 801.15, a glossary of symbols which appear without accompanying text within the product labeling is provided below.

Symbol	Ref	Title of Symbol	Meaning of Symbol
	5.1.1	Manufacturer	Indicates the medical device manufacturer. ^a
	5.1.2	Authorized representative in the European Community	Indicates the authorized representative in the European Community. ^a
	5.1.3	Date of manufacture	Indicates the date when the medical device was manufactured. ^a
	5.1.4	Use-by date	Indicates the date after which the medical device is not to be used. ^a
	5.1.5	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified. ^a
	5.1.6	Catalog number	Indicates the manufacturer's catalog number so that the medical device can be identified. ^a
	5.1.8	Importer	Indicates the entity importing the medical device into the locale. ^a
	5.2.3	Sterilized using ethylene oxide	Indicates the medical device has been sterilized using ethylene oxide. ^a
	5.2.6	Do not re-sterilize	Indicates a medical device that is not to be re-sterilized. ^a
	5.2.8	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened. ^a
	5.2.11	Single sterile barrier system	Indicates a single sterile barrier system. ^a
	5.3.2	Keep away from sunlight	Indicates a medical device that needs protection from light sources. ^a
	5.3.4	Keep dry	Indicates a medical device that needs to be protected from moisture. ^a
	5.3.6	Upper limit of temperature	Indicates the upper limit of temperature to which the medical device can be safely exposed. ^a
	5.4.2	Do not re-use	Indicates a medical device that is intended for one single use only, or for use on a single patient during a single procedure. ^a
	5.4.3	Consult instructions for use www.angiodynamics.com	Indicates the need for the user to consult the instruction for use. ^a
 Cobalt	5.4.10	Contains hazardous substances	Indicates a medical device that contains substances that can be carcinogenic, mutagenic, reprotoxic (CMR) or substances with endocrine disrupting properties. ^a Contains cobalt as a component of stainless steel. This device is not intended for use in the stomach. Exposure of the stainless steel to highly acidic fluids such as gastric fluid can result in leaching of the cobalt from the stainless steel. Cobalt is listed in EC 1272/2008 as a carcinogen class 1B and a reproductive toxin class 1B.
	5.7.7	Medical device	Indicates the items is a medical device. ^a
	5.7.10	Unique Device Identifier	Indicates a carrier that contains unique device identifier information. ^a
	NA	Rx only	Caution: (US) Federal law restricts this device to sale by or on the order of a licensed practitioner. ^b
	NA	Universal Product Number	A Universal Product Number (UPN) code represents the manufacturer's number for an item.
	NA	Quantity in package	To indicate that the adjacent number reflects the number of units contained in the package.
	NA	CE Mark	Manufacturer's declaration of conformity to the Medical Device Regulation 2017/745. ^c

Symbol	Ref	Title of Symbol	Meaning of Symbol
	1135	Recyclable Package	Recyclable Package. ^d
<p>a. EN ISO 15223-1 - Medical Devices – Symbols to be used with medical device labels, labeling and information to be supplied.</p> <p>b. 21 CFR 801.109 - Code of Federal Regulations.</p> <p>c. EU 2017/745 Medical Device Regulations published 5 May 2017</p> <p>d. EN ISO 14021 Environmental labels and declarations. Self-declared environmental claims (Type II environmental labeling)</p>			

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AngioDynamics UK Ltd.