SmartPort⁺ SmartPort Plastic

| | |
|---|-------------------------------|
|  | Instructions For Use 1 |
|  | Instrucciones de uso 22 |
|  | Οδηγίες χρήσης 46 |
|  | Instrukcja używania 68 |
|  | Istruzioni per l'uso 89 |

TABLE OF CONTENTS

| | |
|---|-----------|
| WARNING | 2 |
| DEVICE DESCRIPTION | 2 |
| PORT MATERIAL INFORMATION | 4 |
| INDICATIONS FOR USE | 5 |
| WARNINGS | 7 |
| PRECAUTIONS | 8 |
| MRI SAFETY INFORMATION | 9 |
| NEEDLES | 9 |
| HYPERBARIC CHAMBER SAFETY INFORMATION | 10 |
| POTENTIAL COMPLICATIONS | 10 |
| OPERATIONAL INSTRUCTIONS | 11 |
| PINCH-OFF PREVENTION | 11 |
| IMPLANTATION PREPARATION | 11 |
| PERCUTANEOUS PROCEDURE | 12 |
| CUT-DOWN PROCEDURE | 12 |
| PEELABLE INTRODUCER INSTRUCTIONS | 13 |
| CATHETER TUNNELING PROCEDURE | 14 |
| Detached Catheters | 14 |
| Attached Catheters | 14 |
| CATHETER TO PORT CONNECTION FOR DETACHED CATHETERS | 14 |
| Connecting the Catheter..... | 15 |
| POSITION PORT AND CLOSE INCISION SITE | 15 |
| POWER INJECTION PROCEDURE | 15 |
| Additional Power Injection Information..... | 17 |
| USE AND MAINTENANCE | 17 |
| DETERMINING PORT SYSTEM VOLUMES FOR PORT LOCK PROCEDURES | 18 |
| Equipment | 18 |
| Procedure | 18 |
| SYSTEM OCCLUSION MANAGEMENT | 19 |
| DISCONTINUING USE | 19 |
| WARRANTY | 19 |
| REFERENCES | 20 |

SmartPort⁺

SmartPort Plastic

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Any serious incident which has occurred with the use of this device should be reported to AngioDynamics at complaints@angiodynamics.com and to the National Competent Authority.

Refer to the following web address for contact information for the Competent Authorities. https://ec.europa.eu/health/sites/health/files/md_sector/docs/md_vigilance_contact_points.pdf

The instructions for use are available electronically at ifu.angiodynamics.com

For a copy of the Summary of Safety and Clinical Performance for these devices, please review Eudamed (ec.europa.eu/tools/eudamed) referencing Basic UDI-DI # 50516840024U or contact AngioDynamics Customer Service at +1 800-772-6446

This implantable device is provided with a patient guide, patient implant card, and implant card instructions. The implanting clinician is responsible for reviewing the patient guide with the patient. The implanting clinician is also required to complete the information on the implant card, and provide the completed implant card and patient guide to the patient. The implanting clinician should discuss the risks associated with the devices with the patient.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your sales representative. Inspect prior to use to verify that no damage has occurred during shipping.

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

The SmartPort⁺ and SmartPort^{*} Plastic implantable ports are to be treated as contaminated biomedical waste subsequent to use. Used or unused devices should be disposed of in accordance with hospital, administrative and/or local government policy for such items. Sharps such as the blunt needle and introducer needle should be disposed of in a sharps container.

Uncontaminated device packaging should be recycled if applicable, or disposed of as common waste in accordance with hospital, administrative and/or local government policy for such items.

DEVICE DESCRIPTION

The SmartPort⁺ device with Endexo and *Vortex*^{*} Technology and the SmartPort Plastic device with Vortex Technology are implantable venous access devices designed for repeated access to the vascular system. The SmartPort⁺ and SmartPort Plastic devices are subcutaneous implant devices with one reservoir. The ports are accessed using a non-coring needle which is passed through the skin and into the self-sealing silicone septum covering the reservoir. When used with power injectable needles, the port can be used for power injection of contrast media and contrast enhanced computed tomography (CECT).

Available in plastic or titanium, the SmartPort⁺ device has several configurations of standard, low profile, and mini port body configurations. The ports are offered with 5F, 6F or 8F single lumen *BioFlo*^{*} catheter with Endexo technology that has been shown to be effective in reducing thrombus accumulation (based on platelet count). Reduction of thrombus accumulation was evaluated using acute in vitro models. Pre-clinical in vitro evaluations do not necessarily predict clinical performance with respect to thrombus formation. Vortex technology allows fluid to reach all surfaces in the chamber, helping eliminate dead spaces, resist sludge build-up, and reduce occlusions. The outlet of the vortex port chamber is set at a tangent rather than perpendicularly, allowing flushing action within the port to hyper cleanse the entire chamber leading to decreased sludge build-up and a reduced rate of occlusions.¹

The following components are included with the SmartPort[®]:

- Single Lumen Titanium or Single Lumen Plastic Port
- Single Lumen Catheter (BioFlo, radiopaque, attached or detached)
- Radiopaque Catheter Locking Collar
- 18G Introducer Needle
- 0.038 in Guidewire (J-tipped)
- Peelable Sheath Introducer (Valved or Non-Valved)
- Blunt Needle
- 22G non-coring Needles (1-straight and 1-90°)
- Malleable Tunneler (metal)
- Vein Pick

The SmartPort Plastic device with Vortex technology is offered in a low-profile port body configuration with a 6F or 8F polyurethane catheter.

The following components are included with the SmartPort Plastic:

- Single Lumen Plastic Port
- Single Lumen Catheter (Polyurethane, radiopaque, attached or detached)
- Radiopaque Catheter Locking Collar
- 18G Introducer Needle
- 0.038 in Guidewire (J-tipped)
- Peelable Sheath Introducer (Valved or Non-Valved)
- Blunt Needle
- 22G non-coring Needles (1-straight and 1-90°)
- Malleable Tunneler (metal)
- Vein Pick

In addition to the components above, the SmartPort[®] and SmartPort Plastic devices are supplied with the following:

- Instructions for Use (IFU)
- Product Stickers
- Patient Information Packet containing the Patient Guide, Patient Implant Card and Reminder Band.

The Patient Guide, completed Patient Implant Card and Reminder Band are to be provided to the patient.

The ports are available with either silicone filled or non-filled suture fixation holes. If desired, the suture fixation holes can be used to anchor the port to the subcutaneous tissue. Ports with silicone filled suture holes are intended to prevent tissue in-growth to the suture holes. If needed, filled suture holes are easily accessed through the silicone. All port configurations have a radiopaque identifier (CT mark). The radiopaque catheter is marked at every centimeter and can be cut to the desired length.

In benchtop testing, the SmartPort[®] and SmartPort Plastic devices demonstrated overall compatibility when exposed for clinically relevant durations at clinically relevant dosages to the following agents: contrast medium; antibiotic agents; analgesics; anticoagulation agents; nuclear medicine; and chemotherapy agents inclusive of anti-metabolites, plant alkaloids, antitumor antibiotics, monoclonal antibodies, and antineoplastic agents.

The following commonly used non-coring Needles have been tested with the SmartPort[®] and SmartPort Plastic devices and found no coring of the port septum:

- BD (formally Bard) Power Loc Max 19G x 1.0 in
- BD (formally Bard) Power Loc Max 20G x 1.0 in
- BD (formally Bard) MiniLoc 19G x 1.0 in
- BD (formally Bard) MiniLoc 20G x 1.0 in
- ICU Medical (Formally Smiths Medical) Gripper Plus Power P.A.C. 19G x 1.0 in
- ICU Medical (Formally Smiths Medical) Gripper Plus Power P.A.C. 20G x 1.0 in
- ICU Medical (Formally Smiths Medical) Gripper Plus 19G x 1.0 in
- ICU Medical (Formally Smiths Medical) Gripper Plus 20G x 1.0 in
- PFM Medical EZ Huber 19G x 1.0 in
- PFM Medical EZ Huber 20G x 1.0 in
- Kawasumi K-Shield 19G x 1.0 in
- Kawasumi K-Shield 20G x 1.0 in
- AngioDynamics LifeGuard[®] 19G x 1.0 in
- AngioDynamics LifeGuard 20G x 0.75 in
- PFM Medical Jetcan 22G Huber Needle (Taken from PFM Xcela Port Kit)
- Spectra Medical 22G Huber Needle

If using other manufacturers' needles, please refer to the FDA website for safety information and recommendations at <https://www.fda.gov/MedicalDevices/Safety/default.htm>.

SmartPort®:

- The port body - a small, hollow chamber made of titanium or plastic with a silicone disk (or septum) at the surface for accessing. The port body is offered in various sizes and shapes - "Mini", "Low Profile" and "Standard". The septum is made to be accessed with a special needle (non-coring needle) that will allow sealing of the septum upon removal.
- The catheter - a long, soft flexible tube made of polyurethane with Endexo Technology. One end of the catheter is securely connected to the port, and the other end is placed into a large vein in your chest.

SmartPort Plastic:

- The port body - made of plastic with a silicone disk (or septum) at the surface for accessing. Offered as "Low Profile" port.
- The catheter - a long, soft flexible tube made of polyurethane.

PORT MATERIAL INFORMATION

| Component | Material | Type of Body Contact | Weight of material / Weight of implanted device (in g) |
|--|---|----------------------------------|--|
| Septum Mini Titanium | Medical grade silicone elastomer, (dimethyl siloxane) | Blood path indirect, Tissue/bone | 0.9/10.6 |
| Septum Low profile Titanium | Medical grade silicone elastomer, (dimethyl siloxane) | | 0.9/13.2 |
| Septum Standard Titanium port | Medical grade silicone elastomer, (dimethyl siloxane) | | 0.9/15.4 |
| Septum Plastic port | Medical grade silicone elastomer, (dimethyl siloxane) | | 0.9/7.7 |
| O-ring gasket (Standard Titanium ports only) | Medical grade silicone elastomer, (dimethyl siloxane) | Blood path indirect, Tissue/bone | Negligible* |
| Adhesive (Plastic ports) | Epoxy/amide polymer | Blood path indirect, Tissue/bone | Negligible* |
| Adhesive (Titanium ports) | Medical grade silicone elastomer, (dimethyl siloxane) | | Negligible* |
| Suture Feet Plugs (all ports with filled fixation holes) | Medical grade silicone elastomer, (dimethyl siloxane) | Blood path indirect, Tissue/bone | Negligible* |
| Catheter tubing ink SmartPort® | Titanium dioxide, and aluminum oxide pigments with a binder of 2-(2-butoxyethoxy)ethyl dihydrogen phosphate, N,Ndimethylcyclohexanamine | Blood path direct, Tissue/bone | Negligible* |
| Catheter tubing ink SmartPort Plastic | Tall oil esters | | Negligible* |
| Stem (Mini Titanium port) | Alloy containing 90% Titanium, 6% Aluminum, 4% Vanadium | Blood path indirect, Tissue/bone | 0.2/10.6 |
| Stem (Low Profile Titanium port) | Alloy containing 90% Titanium, 6% Aluminum, 4% Vanadium | | 0.2/13.2 |
| Stem (Standard Titanium port) | Alloy containing 90% Titanium, 6% Aluminum, 4% Vanadium | | 0.2/15.4 |
| Stem (Low Profile Plastic port) | Alloy containing 90% Titanium, 6% Aluminum, 4% Vanadium | | 0.2/7.7 |

| Component | Material | Type of Body Contact | Weight of material / Weight of implanted device (in g) |
|--|---|----------------------------------|--|
| Catheter shaft in SmartPort [®] Mini Titanium port | 20% Aliphatic polyurethane resin, 8% barium sulfate, 0.6% fluorinated poly-ester-urethane, 0.06% copper phthalocyanine pigment | Blood path direct, Tissue/bone | 3.1/10.6 |
| Catheter shaft in SmartPort [®] Low Profile Titanium port | 16% Aliphatic polyurethane resin, 7% barium sulfate, 0.5% fluorinated poly-ester-urethane, 0.05% copper phthalocyanine pigment | | 3.1/13.2 |
| Catheter shaft in SmartPort [®] Standard Titanium port | 14% Aliphatic polyurethane resin, 6% barium sulfate, 0.4% fluorinated poly-ester-urethane, 0.04% copper phthalocyanine pigment | | 3.1/15.4 |
| Catheter shaft in SmartPort [®] Low Profile Plastic port | 27% Aliphatic polyurethane resin, 12% barium sulfate, 0.8% fluorinated poly-ester-urethane, 0.08% copper phthalocyanine pigment | | 3.1/7.7 |
| Catheter shaft in SmartPort Plastic Low Profile Plastic port | 28% Aliphatic polyurethane resin, 12% barium sulfate | | 3.1/7.7 |
| Plastic collar (Standard Titanium port) | 0.6% polysulfone resin, 0.09% barium sulfate | Blood path indirect, Tissue/bone | 0.1/15.4 |
| Plastic collar (Low Profile Titanium port) | 0.7% polysulfone resin, 0.1% barium sulfate | | 0.1/13.2 |
| Plastic collar (Low Profile Plastic ports) | 1.2% polysulfone resin, 0.18% barium sulfate | | 0.1/7.7 |
| Titanium Collar (Mini Titanium ports) | Alloy containing 90% Titanium, 6% Aluminum, 4% Vanadium | | 0.3/10.6 |
| Port housing Standard Titanium port | Alloy containing 90% Titanium, 6% Aluminum, 4% Vanadium | Blood path indirect, Tissue/bone | 11.1/15.4 |
| Port housing Mini Titanium port | Alloy containing 90% Titanium, 6% Aluminum, 4% Vanadium | | 6.1/10.6 |
| Port housing Low Profile Titanium port | Alloy containing 90% Titanium, 6% Aluminum, 4% Vanadium | | 8.8/13.2 |
| Port housing Low Profile Plastic port | 40% polysulfone resin, 3% tungsten, 1% [2-(2-butoxyethoxy) ethyl dihydrogen phosphate, N,Ndimethylcyclohexanamine] | | 3.4/7.7 |

* NOTE: "Negligible" means the weight of this component does not appreciably affect overall weight of the device. Component is less than 0.1 g in weight.

INDICATIONS FOR USE

The ports are indicated for patients who require long-term access to the central venous system for blood specimen withdrawal and administration of fluids including but not limited to hydration fluids, chemotherapy, analgesics, nutritional therapy and blood products, as well as the administration and adequate removal of nuclear medicine.

When used with power injectable needles, the ports are indicated for power injection of contrast media. For power injection of contrast media, the maximum recommended infusion rate is 5 mL/s with 19G or 20G non-coring power injectable needles or 2 mL/s with a 22G non-coring power injectable needle.

| Needle Size (G), non-coring power injectable | Catheter Size (F) | Maximum Recommended Flow Rate Setting (mL/s) | Maximum Recommended Pressure Setting (psi) |
|--|-------------------|---|---|
| 19/20 | 5, 6, and 8 | 5 | 300 |
| 22 | 5, 6, and 8 | 2 | 300 |

Considering the indication for long-term access to the vascular system, the ports may, at the discretion of the clinician, be employed in the clinical management, evaluation, and monitoring of conditions such as, but not limited to, the following:

- Dehydration
- Malnutrition
- Malignant conditions
- Fluid and electrolyte imbalances
- Pain syndromes
- Anemia
- Hematologic abnormalities

INTENDED USE:

The SmartPort[®] and SmartPort Plastic implantable ports systems are intended for subcutaneous implantation for long-term vascular access. When used with power injectable needles, the port is indicated for power injection of contrast media.

INTENDED PATIENT POPULATION:

These ports can be used in adult patients who require long-term vascular access for diagnostic and therapeutic management of their disease states. Clinicians must consider variations in individual patient anatomy and physiology due to size and age. When determining patient selection, catheter diameter and port body size should be considered. Appropriate guidance, vein assessment and insertion techniques for implanted port placement should be employed.

INTENDED USERS:

The ports should be inserted, manipulated, and removed by a qualified, licensed implanting clinician or other qualified health care professional under the direction of an implanting clinician.

These ports should be used only by implanting clinicians and healthcare professionals already trained in port access, care and maintenance. Users may receive additional product training, as needed, from a clinical team representative from AngioDynamics by contacting customer service at +1800-772-6446.

INTENDED CLINICAL BENEFITS:

- SmartPort[®] and SmartPort Plastic implantable ports offer an option to meet long-term vascular access needs.
- The ports are indicated for long-term use (> 30 days^{2,3}) and are to be removed as clinically indicated⁴.
- SmartPort[®] and SmartPort Plastic implantable ports offer the potential benefit of low risk for device-associated infection⁴.

CONTRAINDICATIONS

- Catheter insertion in the subclavian vein medial to the border of the first rib, an area which is associated with higher rates for pinch-off.⁵
- Presence of infection, bacteremia, or septicemia.
- Past irradiation of prospective insertion site.
- Previous episodes of venous thrombosis or vascular surgical procedures at the prospective placement site.
- Local tissue factors to prevent proper device stabilization and/or access.
- Hypercoagulopathy unless considerations are made to place the patient on anticoagulation therapy.
- Presence or suspicion of allergic reaction to materials contained in this device.
- Anatomy is insufficient to accommodate size of the port or the catheter.
- Demonstrated intolerance for an implanted device.

WARNINGS

- The device is to be implanted, used, maintained and removed in accordance with current evidence-based guidance for infection prevention from agencies/organizations such as Centers of Disease Control (CDC) or the World Health Organization. It is recommended that institutional policies and procedures are updated periodically to reflect current evidence-based guidelines for infection prevention.
- During placement through a non-valved introducer sheath, hold thumb over the exposed opening of sheath to prevent air embolism or patient injury may occur. The risk of air embolism is reduced by performing this part of the procedure with the patient performing the Valsalva maneuver.
- Do not suture catheter to port, port stem, or surrounding tissue. Any damage or constriction of catheter may compromise power injection performance and catheter integrity.
- Do not use syringes smaller than 10 mL syringe when accessing the port as system damage can occur. Flushing occluded catheters with small syringes can create excessive pressures within the port system.
- Do not forcefully flush the port system with any syringe size. After confirmation of patency by detecting no resistance and the presence of a blood return, use syringes appropriately sized for the medication being injected. Do not transfer the medication to a larger syringe.⁵
- Failure to use a power injectable needle with the SmartPort® or SmartPort Plastic device for a power injection procedure may result in port system failure and patient injury may occur. Refer to the Power Injection Procedure section in this booklet for additional information and instructions.
- Failure to ensure patency of the catheter prior to power injection studies may result in port system failure and patient injury may occur.
- Do not power inject through a port system that exhibits signs of clavicle-first rib compression or pinch-off as it may result in port system failure and patient injury may occur.
- Failure to warm contrast media to body temperature may result in port system failure and patient injury may occur.
- Do not exceed 300 psi (2068 kPa) pressure limit setting or the maximum recommended flow rate setting. Exceeding the maximum flow rate may result in port system failure and/or catheter tip displacement and patient injury may occur.
- If local pain, swelling or signs of extravasation are noted, the injection should be stopped immediately as patient injury may occur.
- Absence of a blood return or a poor blood return can be a sign of a potential complication such as occlusion, kinking, breakage, Pinch-Off Syndrome, fibrin formation, thrombosis or malposition. This should be evaluated prior to device usage.
- A blood return should be present prior to usage of device for any therapy or testing.
- Do not attempt to measure the patient's blood pressure on the arm in which a peripheral system is located, since catheter occlusion or other damage to the catheter could occur.
- If the patient complains of pain, or there is swelling when the device is flushed or when medication or contrast media is administered, evaluate the device for infiltration, proper needle placement, and potential complications such as occlusion, kinking, breakage, Pinch-Off Syndrome, thrombosis or malposition. Failure to assess these complaints or observations can lead to device failure and patient injury may occur.
- Power injection machine pressure limiting (safety cut-off) settings may not prevent over-pressurization of an occluded device.
- Reuse of single-use devices creates a potential risk of patient or user infections.
- Contamination of the device may lead to injury, illness or death of the patient.
- Reprocessing may compromise the integrity of the device and/or lead to device failure and patient injury may occur.
- Review all fluids to be used with this product for any incompatibility with plastics or plastic adhesives.
- 18 Gauge Introducer needle, Non-coring needle, Guidewire, Tunneler contain cobalt. Cobalt is classified as CMR 1B and is present in a concentration above 0.1% weight by weight.

PRECAUTIONS

- Carefully read and follow all instructions prior to use.
- Only licensed health care practitioners should insert, manipulate, and remove these devices.
- Prime the device with sterile heparinized saline or 0.9% Sodium Chloride for Injection solution to help avoid air embolism.
- Some patients may be hypersensitive to heparin or suffer from heparin induced thrombocytopenia (HIT) and these patients must not have their port primed with heparinized saline.
- When using an introducer kit, verify that the catheter fits through the introducer sheath.
- When using percutaneous introducers:
 - Carefully insert the introducer over the guidewire to avoid inadvertent penetration to vital structures in the thorax.
 - To avoid blood vessel damage, do not allow the percutaneous introducer sheath to remain indwelling in the blood vessel without the internal support of a catheter or dilator.
 - Simultaneously advance the sheath and dilator with rotational motion to help prevent sheath damage.
- To avert device damage and/or patient injury during catheter placement:
 - Avoid accidental device contact with sharp instruments and mechanical damage to the catheter material.
 - Use only smooth-edged atraumatic clamps or forceps.
 - Do not use the catheter if there is any evidence of mechanical damage or leaking.
 - Avoid sharp or acute angles during implantation which could compromise the patency of the catheter lumen.
 - Carefully follow the Catheter to Port Connection section of this DFU for instructions to ensure proper device connection and to avoid catheter damage.
- Assure tight connection between port body and catheter.
- After implantation or usage of the port, the system should be flushed to clear infusates and/or blood components in order to maintain patency, per institutional protocol. Refer to the Use and Maintenance section of this DFU for recommendations.
- For precise drug administration please refer to the individual pharmaceutical instructions.
- Prior to any treatment palpate correct position of the port body and assure no signs or symptoms of port site irritation or infection exist.
- Only use non-coring needles to access the port septum. The non-coring needle tip is integral to prevent port septum damage.
- Palpate the port and port septum then access the septum with a non-coring needle at a 90° angle.
- Puncture skin directly over septum and gently advance needle through septum until it contacts bottom of port chamber. Do not apply excessive force once the needle contacts the port floor.
- Prior to injection or infusion, aspirate to ensure a brisk blood return. If blood return is not realized, see System Occlusion Management section below.
- If more than one drug is to be administered, between the individual drug applications, flush the system with 5 to 10 mL normal saline for injection to prevent drug interaction.
- After any infusion, injection or bolus application, the system should be flushed with normal saline. The system should be locked with a heparin solution or per institutional protocol to prevent thrombotic occlusion of the catheter. In all cases, institutional protocols and current clinical practice guidelines and standards for flushing and locking implanted ports should be considered.
- Do not use alcohol to soak or de clot polyurethane catheters because alcohol is known to degrade the polyurethane material over time with repeated and prolonged exposure.

MRI SAFETY INFORMATION

Magnetic Resonance
Conditional

MRI Conditional Information

Non-clinical testing has demonstrated that the SmartPort⁺ and SmartPort Plastic devices are MR Conditional. A patient with these devices can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5-Tesla (1.5 T), 3-Tesla (3 T), 7-Tesla (7 T).
- Maximum spatial field gradient of 4,500 G/cm (45 T/m).
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg (Normal Operating Mode) for 15 minutes of scanning.
- At 7-T the device must remain outside the transmit coil.

RF heating Under the scan conditions defined above, the SmartPort⁺ and SmartPort Plastic devices are expected to produce a maximum temperature rise of less than or equal to 4 °C after 15 minutes of continuous scanning.

MR Artifact

In non-clinical testing, the image artifact caused by the device extends approximately 26 mm from the SmartPort⁺ titanium configurations and 13 mm from the SmartPort⁺ plastic and SmartPort Plastic configurations when imaged with a spin-echo or gradient-echo pulse sequence in a 3-T MRI system.

IMPORTANT NOTE: The 7-Tesla MR system (MAGNETOM, Terra, Siemens Medical Solutions USA, Inc.) currently has a transmit/receive RF head and a transmit/receive RF knee coil. There is no transmit/receive RF body coil. During the intended use of the SmartPort⁺ and SmartPort Plastic, the Device is not implanted in the anatomical regions of the head or the knee. Therefore, the MRI testing performed on this product was limited to an assessment of magnetic field interactions, only. Once a transmit/receive RF body coil becomes commercially available for this clinical scanner, it will be necessary to conduct further MRI testing to evaluate MRI-related heating and artifacts for the SmartPort⁺ and SmartPort Plastic devices.

Precaution: The RF heating behavior does not scale with static field strength. Devices which do not exhibit detectable heating at one field strength may exhibit high values of localized heating at another field strength.

NEEDLES

Use of non-coring needles (19G, 20G or 22G) is recommended for all procedures. Septum puncture life for SmartPort⁺ and SmartPort Plastic devices includes:

| Needle Size (G), non-coring power injectable | Maximum Needle Sticks |
|--|-----------------------|
| 19/20 | 1000 |
| 22 | 1500 |

HYPERBARIC CHAMBER SAFETY INFORMATION

The ports were exposed to 3 ATA (Atmosphere Absolute) in a hyperbaric chamber followed by an emergency blow procedure and no gas escape or surface tension changes were observed. Benchtop testing was performed following hyperbaric chamber exposure to ensure product functionality (see burst results below*).

| Port Configuration | Catheter Size | Catheter Material | Max Hyperbaric Chamber Pressure | Average Static Burst Pressure PSI (kPa) ^{††} | Static Burst Pressure Range PSI (kPa) ^{††} |
|---------------------|----------------|-------------------|---------------------------------|---|---|
| Standard Ti | 8F (2.6 mm) | BioFlo | 3 ATA | 246 (1695) | 220 - 265 (1517 - 1827) |
| Mini Ti | 5F (1.8 mm) | BioFlo | 3 ATA | 218 (1506) | 167 – 316 (1151 - 2179) |
| Mini Ti | 6F (2.1 mm) | BioFlo | 3 ATA | 282 (1943) | 261 - 298 (1800 - 2055) |
| Low Profile Plastic | 6F (2.1 mm) | Std. PU | 3 ATA | 184 (1266) | 157 - 202 (1083 - 1393) |
| Low Profile Ti | 8F (2.6 mm) | Std. PU | 3 ATA | 279 (1921) | 264 - 298 (1820 - 2055) |

* Benchtop testing, devices were not tested clinically for performance after hyperbaric chamber exposure.

†† Worst case static burst pressure of port-catheter assembly specified.

POTENTIAL COMPLICATIONS

- Air or catheter (or catheter fragments) embolism
- Allergic reaction
- Arterial puncture
- Arteriovenous fistula
- Bacteremia
- Bleeding
- Brachial plexus injury
- Cardiac arrhythmia
- Cardiac puncture
- Cardiac tamponade
- Catheter or port erosion through skin/vessel
- Catheter occlusion, malposition, dislodgment, fragmentation, migration, disconnection or rupture
- Catheter occlusion or breakage caused by pinching between clavicle and first rib
- Catheter thrombosis
- Chylothorax
- Death
- Drug extravasation
- Endocarditis
- Fibrin sheath formation
- Hematoma
- Hemothorax
- Hydrothorax
- Inflammation
- Infection
- Implant rejection
- Implant rotation or extrusion
- Implantation site necrosis or scarring of skin over implant area
- Laceration or perforation of vessel
- Necrosis or scarring over skin implant area
- Nerve damage
- Pain at or around port pocket site
- Peritonitis
- Pneumothorax
- Risks normally associated with local and general anesthesia, surgery, and post-operative recovery
- Stuck or difficult to remove catheters
- Thoracic duct injury
- Thromboembolism
- Thrombophlebitis
- Tunnel infection
- Twiddler's syndrome
- Vascular thrombosis
- Vessel trauma

HOW SUPPLIED

Contents supplied STERILE using a Ethylene Oxide (EO) process. Store in a cool, dry, dark place. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

OPERATIONAL INSTRUCTIONS

Clinical signs of catheter pinch-off may include, but not limited to:

- difficulty with blood withdrawal
- resistance to infusion of fluids
- patient position changes required for infusion of fluids or blood withdrawal


Radiologic signs of catheter pinch-off:

- Grade 1 or 2 distortion on chest X-ray. Pinch-off should be evaluated for degree of severity prior to explantation. Patients indicating any degree of catheter distortion at the clavicle/1st rib area should be followed diligently. There are grades of pinch-off that should be recognized with appropriate chest x-ray as follows:⁵

| Grade | Severity | Recommended Action |
|---------|--|---|
| Grade 0 | No Distortion | No Action |
| Grade 1 | Distortion present without luminal narrowing | Chest x-ray should be taken every one (1) to three (3) months to monitor progression of pinch off to grade 2 distortion. Shoulder positioning during chest x-rays should be noted as it can contribute to changes in distortion grades. |
| Grade 2 | Distortion present with luminal narrowing | Removal of catheter should be considered |
| Grade 3 | Catheter transection or fracture | Prompt removal of catheter |

PINCH-OFF PREVENTION

The risk of pinch-off syndrome can be avoided by inserting the catheter via the internal jugular vein (IJ). Catheters placed percutaneously or through a cut-down into the subclavian vein should be inserted at the junction of the outer and middle thirds of the clavicle, lateral to the thoracic outlet. The catheter should not be inserted into the subclavian vein medially, because such placement can lead to compression of the catheter between the first rib and the clavicle, which can cause damage and even catheter fracture. A radiographic confirmation of catheter placement should be made to ensure that the catheter is not being pinched by the first rib and clavicle.

IMPLANTATION PREPARATION

Implantation and catheter insertion may be accomplished by a number of techniques. Choice of the proper procedure is the responsibility of the implanting clinician.

1. Select implantation procedure to be used.

Note: Recommended veins for chest placement are the internal jugular or lateral subclavian. Refer to Pinch-Off Prevention section if inserting the catheter via subclavian vein.

2. Select site for port placement.

Note: Port pocket site selection should allow for port placement in an anatomic area that provides good port stability, does not interfere with patient mobility, create pressure points, or interfere with clothing. Consider the amount of cutaneous tissue over the port septum as excessive tissue will make access difficult. Conversely, too thin a tissue layer may lead to port erosion. A tissue thickness of 0.5 cm to 2 cm is appropriate.

3. Complete patient implant record, documenting device specific information according to institutional policies and procedures including reorder number and lot number.

4. Create sterile field and open tray.
5. Surgically prep and drape the operative site.
6. Perform adequate anesthesia.
7. Using the non-coring needle provided with the port kit, flush catheter and port with 5 mL of 0.9% Sodium Chloride for Injection and Clamp the catheter closed at the section that will be trimmed off prior to attachment to the port body and/or prior to advancing the catheter into the vasculature. The clamped section should be trimmed off before final port placement. For port models with detached catheters, prime the catheter with 0.9% Sodium Chloride for Injection using the blunt needle provided in kit and the port body with the non-coring needle provided.

Note: Clamped catheter segments should be trimmed off before final port placement.

Precautions: Some patients may be hypersensitive to heparin or suffer from heparin induced thrombocytopenia (HIT) and these patients must not have their port primed with sterile heparinized saline. Avoid handling the catheter with sharp objects. When handling the catheter, padded hemostats, vascular clamps, or tube occluding forceps should be used. Instruments with teeth should never be used to grasp the catheter. Damaging the catheter prior to or during insertion may cause catheter fracture in the vessel. The catheter should only be grasped at the end that will be trimmed prior to insertion.

8. Place patient in the Trendelenburg position with head turned away from the intended venipuncture site.

PERCUTANEOUS PROCEDURE


Note: If using a vascular access system other than that included with the SmartPort® or SmartPort Plastic device consult that manufacturer's Instructions for Use.

- 1 Select appropriate vessel for introduction of the catheter.
- 2 Gain access to the selected vein with the introducer needle attached to a syringe.
- 3 Confirm proper placement of the needle tip in the vessel by aspiration of blood and/or by means of ultrasound guidance.
- 4 Remove syringe, leaving needle in vessel.

Warning: Place a finger over the hub of the needle to minimize blood loss and risk of air embolism.

- 5 Straighten "J" tip of guidewire with tip straightener and insert tapered end of tip straightener into the needle. Pass J-tip guidewire through the needle and advance to the superior vena cava. Advance the guidewire as far as appropriate for the procedure.
- 6 Verify correct positioning using fluoroscopy or other appropriate technology.
- 7 Gently withdraw and remove the needle and secure guidewire.

Precaution: If the guidewire must be withdrawn while the needle is inserted, remove both the needle and the wire as a unit to help prevent the needle from damaging or shearing the guidewire.

CUT-DOWN PROCEDURE

1. Use a cut-down incision to expose the entry vein of choice.
2. Perform vessel incision after the vessel is isolated and stabilized to prevent bleeding and air embolism.
3. With the aid of a vein pick, insert the catheter through a small venotomy into the isolated vein and advance catheter tip to the distal third of the superior vena cava, near the cavo-atrial junction or alternate tip location, as clinically indicated.

- Secure placement and withdraw vein pick.
- Verify correct catheter tip position using fluoroscopy or appropriate technology.


Warning: Do not suture catheter to port, port stem, or surrounding tissue. Any damage or constriction of catheter may compromise power injection performance and catheter integrity.

PEELABLE INTRODUCER INSTRUCTIONS

Note: If using an introducer system other than that included with the SmartPort® or SmartPort Plastic device, consult that manufacturer's Instructions for Use.

Precaution: When using an introducer kit, verify that the catheter fits easily through the introducer sheath.

- Flush introducer with saline prior to use.
- Insert the dilator through the sheath and lock in place.
- Make a small incision to facilitate advancement of dilator/sheath assembly over the guidewire, taking care not to allow contact of scalpel with wire.
- Thread the dilator/sheath assembly over the guidewire.

Precaution: Carefully insert the introducer over the guidewire to avoid inadvertent penetration to vital structures in the thorax. To avoid blood vessel damage, do not allow the percutaneous introducer sheath to remain indwelling in the blood vessel without the internal support of a catheter or dilator. Simultaneously advance the sheath and dilator with rotational motion to help prevent sheath damage.

Note: Fluoroscopic observation may be advisable. Attaching a clamp or hemostat to the proximal end of the guidewire will prevent inadvertently advancing the guidewire entirely into the patient.

- Unlock the dilator from the sheath handle by turning the dilator collar counterclockwise.
- Gently remove the dilator and guidewire, leaving the sheath as a conduit to the vessel.


Warning: If using a non-valved introducer sheath, hold thumb over the exposed opening of sheath to prevent air embolism. The risk of air embolism is reduced by performing this part of the procedure with the patient performing the Valsalva maneuver.

- Advance the catheter through the sheath and into the vessel.

Note: To prevent kinking the catheter, it may be necessary to advance in small steps by grasping the catheter close to the sheath. Some resistance may be felt when advancing the catheter.

Precautions: Avoid accidental device contact with sharp instruments and mechanical damage to the catheter material. Use only smooth-edged atraumatic clamps or forceps. Avoid sharp or acute angles during implantation which could compromise the patency of the catheter lumen.

- When the catheter is properly positioned, crack and peel the sheath handle in half and continue to pull so that the sheath separates longitudinally while withdrawing the sheath from the vein. Make sure the catheter is not dislodged from the vessel.

CATHETER TUNNELING PROCEDURE

1. Create a subcutaneous pocket for the port using blunt dissection.

Note: Do a trial placement to verify that the pocket is large enough to accommodate the port and that the port does not lie beneath the incision.

Detached Catheters

1. Insert the tip of the tunneler into a small incision at the venous entry site.
2. Remove the catheter lock from the catheter.
3. Remove the tunneler cap, then attach the tip of the catheter onto the tunneler barbs. Barb threads must be completely covered by the catheter to adequately secure the catheter as it is pulled through the tunnel. A suture may be tied around the catheter between the tunneler body and the large barb to hold it more securely.

Precaution: Avoid inadvertent puncture of the skin or fascia with the tip of the tunneler.

4. Using the tapered tip of the tunneler to perform blunt dissection, create a subcutaneous tunnel from the incision at the venous entry site through to the port pocket site.
5. Gently pull the tunneler through the insertion site to advance catheter into the tunnel.

Note: If resistance is encountered, further blunt dissection may facilitate insertion.

6. When the catheter tip has fully emerged from the tunnel cut the catheter off of the tunneler at approximately 45-degree angle to facilitate placement of locking collar on the catheter.
7. Do not pull the catheter to detach the tunneler as this may cause damage to the catheter.

Attached Catheters

1. Remove the tunneler cap and connect the catheter tip onto the barbed end of tunneler.

Precaution: Avoid inadvertent puncture of the skin or fascia with the tip of the tunneler.

2. Using the tapered tip of the tunneler to perform blunt dissection, create a subcutaneous tunnel starting at the pocket site to the incision at the venous entry site.
3. Gently pull the tunneler through the pocket site to advance catheter into the tunnel.

Note: If resistance is encountered, further blunt dissection may facilitate insertion.


4. When the distal tip has fully emerged from the tunnel to the incision at the venous entry site, trim off the end of the catheter attached to the tunneler. Do not pull the catheter to detach the tunneler as this may cause damage to the catheter.
5. Estimate the catheter length required for the tip placement by placing the catheter on the chest along the venous path to the lower third of the SVC at or near the cavo-atrial junction.
6. Trim the catheter to proper length at a 90° angle allowing sufficient slack to permit body movement, port connection and verify that the catheter is not kinked.

CATHETER TO PORT CONNECTION FOR DETACHED CATHETERS

1. Prime all port components per the instructions in the Implantation Preparation section.
2. Flush all air from the port by using a 10mL or larger syringe with non-coring needle filled with 0.9% Sodium Chloride for Injection. Insert the needle through the septum and inject the fluid while pointing the stem up.

Connecting the Catheter

- Place catheter lock onto catheter. The catheter lock is bi-directional and can be mounted on the catheter in either direction.
- Trim the catheter to proper length at a 90° angle allowing sufficient slack to permit body movement, port connection and verify that the catheter is not kinked.
- Advance catheter over port stem to the midway point, slightly past the barb.
- Advance catheter lock until it engages with tactile and/or audible feedback.

Note: To ensure proper assembly of port and catheter lock connector, a minimal gap (less than 0.5 mm) is expected.

Note: If the catheter and locking collar are connected and then disconnected, the catheter proximal end must be re-cut to ensure a safe re-connection to the port.

Precaution: Prior to advancing the catheter lock connector, ensure that the catheter is properly positioned. A catheter not advanced to the proper region may not seat securely and lead to dislodgement and extravasation. The catheter must be straight with no sign of kinking. A slight pull on the catheter is sufficient to straighten it. Advancing the catheter lock over a kinked catheter may damage the catheter.

POSITION PORT AND CLOSE INCISION SITE

- Place the port in the subcutaneous pocket away from the incision line.
- Verify correct tip position using fluoroscopy or other technology.
- Access the port with a non-coring needle and assess patency. Conduct flow studies on the catheter using a non-coring needle and a syringe to confirm that the flow is not obstructed, that no leak exists, and that the catheter is correctly positioned.
- Aspirate to confirm the ability to draw blood.
- Secure to the underlying fascia using one non-absorbable, monofilament suture per suture hole.
- Close the incision site(s).
- Access with non-coring needle to confirm patency by aspirating blood and flushing.
- After implantation or usage of the port, the system should be flushed to clear infusates and/or blood components and locked in order to maintain patency. Refer to the Use and Maintenance section of this DFU for recommendations.

Precaution: Some patients may be hypersensitive to heparin or suffer from heparin induced thrombocytopenia (HIT) and these patients must not have their port primed with heparinized saline.

POWER INJECTION PROCEDURE

- Verify that the patient has a SmartPort⁺ or SmartPort Plastic device by of the following means:
 - Check patient's chart/implant record for a SmartPort⁺ or SmartPort Plastic device or product sticker.
 - SmartPort⁺ and SmartPort Plastic power-injectable ports can be identified by the Smart Angle[®] technology on the CT and CT Low-Profile models. The CT engraving on all models may be identified through fluoroscopy, chest x-ray, or CT scout scan. Each SmartPort⁺ or SmartPort Plastic kit contains an education packet including an information booklet, implant card, implant card leaflet, wrist band, and patient guide that should be provided to the patient or caregiver following port implantation.

Precaution: Prior to any treatment palpate correct position of the port body and assure no signs or symptoms of port site irritation or infection exist.

Precaution: Only use non-coring needles to access the port septum. The special needle tip prevents damage of the port septum.

Precaution: Palpate the port and port septum, then access the septum with a non-coring needle at a 90-degree angle.

Precaution: Puncture skin directly over septum and gently advance needle through septum until it contacts bottom of port chamber. Do not apply excessive force once the needle contacts the bottom of port chamber.

2. Access the SmartPort® or SmartPort Plastic device with a power injectable needle.

Warning: Failure to use a power injectable needle with the SmartPort® or SmartPort Plastic device for a power injection procedure may result in port system failure and patient injury may occur.

Note: Consult manufacturer's Instructions for Use for power injectable needle.

3. Attach a syringe filled with sterile normal saline.
4. Check for patency with the patient in the position that they will assume during the CECT procedure. If possible, the patient should receive power injection with the patient's arm vertically above the shoulder with the palm of the hand on the face of the gantry during injection. This allows for uninterrupted passage of injected contrast through the axillary and subclavian veins at the thoracic outlet.
5. Aspirate for adequate blood return and vigorously flush the port with at least 10 mL of 0.9% Sodium Chloride for Injection.

Warning: Failure to ensure patency of the catheter prior to power injection studies may result in port system failure and patient injury may occur. Do not power inject through a port system that exhibits signs of clavicle-first rib compression or pinch-off as it may result in port system failure and patient injury may occur.

6. Detach syringe.
7. Warm contrast media to body temperature.

Warning: Failure to warm contrast media to body temperature prior to power injection may result in port system failure and patient injury may occur.

Note: Follow institutional protocol to verify correct catheter tip position prior to power injection.

8. Attach the power injection device to the power injectable needle ensuring connection is secured.
9. Check table below to confirm the maximum flow rate and maximum pressure setting.

| Power Injectable Non-Coring Needle Size | 19G | 20G | 22G |
|---|--------------------|--------------------|--------------------|
| Maximum Flow Rate Setting | 5 mL/sec | 5 mL/sec | 2 mL/sec |
| Maximum Pressure Setting | 300 psi (2068 kPa) | 300 psi (2068 kPa) | 300 psi (2068 kPa) |

Warning: Do not exceed 300 psi (2068 kPa) pressure limit setting or the maximum recommended flow rate setting. Exceeding the maximum flow rate may result in port system failure and/or catheter tip displacement, and patient injury may occur.

10. Instruct the patient to communicate immediately any pain or change in feeling during the injection.
11. Inject warmed contrast, taking care not to exceed the flow rate limits.

Warning: If local pain, swelling or signs of extravasation are noted, the injection should be stopped immediately as patient injury may occur.

12. Disconnect the power injection device.

13. After usage of the port, the system should be flushed to clear infusates and/or blood components and locked in order to maintain patency. Refer to the Use and Maintenance section of this DFU for recommendations.

Precaution: Some patients may be hypersensitive to heparin or suffer from heparin induced thrombocytopenia (HIT) and these patients must not have their port primed with heparinized saline.

Additional Power Injection Information

| Port Configuration | Catheter Size | Catheter Material | Max. Flow Rate Limit Setting | Average Flow Rate Performance [†] | Average Static Burst Pressure PSI (kPa) ^{††} | Static Burst Pressure Range PSI (kPa) ^{††} |
|---------------------|---------------|-------------------|------------------------------|--|---|---|
| Standard Ti | 5F (1.8 mm) | BioFlo | 5 mL/s | 5 mL/s | 221 (1521) | 195 - 232 (1344 - 1600) |
| Standard Ti | 6F (2.1 mm) | BioFlo | 5 mL/s | 5 mL/s | 210 (1448) | 193 - 231 (1331 - 1593) |
| Standard Ti | 8F (2.6 mm) | BioFlo | 5 mL/s | 5 mL/s | 224 (1547) | 178 - 253 (1227 - 1744) |
| Low Profile Ti | 5F (1.8 mm) | BioFlo | 5 mL/s | 5 mL/s | 281 (1935) | 219 - 289 (1510 - 1993) |
| Low Profile Ti | 8F (2.6 mm) | BioFlo | 5 mL/s | 5 mL/s | 242 (1671) | 192 - 262 (1324 - 1806) |
| Mini Ti | 5F (1.8 mm) | BioFlo | 5 mL/s | 5 mL/s | 271 (1871) | 261 - 281 (1800 - 1937) |
| Mini Ti | 6F (2.1 mm) | BioFlo | 5 mL/s | 5 mL/s | 264 (1819) | 246 - 279 (1696 - 1924) |
| Low Profile Plastic | 5F (1.8 mm) | BioFlo | 5 mL/s | 5 mL/s | 169 (1165) | 162 - 176 (1117 - 1213) |
| Low Profile Plastic | 8F (2.6 mm) | BioFlo | 5 mL/s | 5 mL/s | 157 (1085) | 137 - 167 (945 - 1151) |
| Low Profile Plastic | 6F (2.1 mm) | Std. PU | 5 mL/s | 5 mL/s | 176 (1216) | 158 - 194 (1089 - 1338) |
| Low Profile Plastic | 8F (2.6 mm) | Std. PU | 5 mL/s | 5 mL/s | 159 (1094) | 150 - 165 (1034 - 1138) |

[†] Represents flow capability of port and catheter assembly for power injection of contrast media with 19G infusion set.

^{††} Worst case static burst pressure of port-catheter assembly specified.

Note: The pressures provided are worst case for the port/catheter configuration specified.

Note: CT injector pressure limit should be set at a maximum of 300 psi (2068 kPa).

USE AND MAINTENANCE

To help prevent clot formation and catheter blockage, the SmartPort[®] or SmartPort Plastic should be flushed and locked to clear infusates and/or blood components in order to maintain patency. If the port remains unused for long periods of time, the lock should be changed at least once every four weeks or per institutional protocol. Refer to the tables below for manufacturer recommendations.

Precaution: Some patients may be hypersensitive to heparin or suffer from heparin induced thrombocytopenia (HIT) and these patients must not have their port primed with heparinized saline.

DETERMINING PORT SYSTEM VOLUMES FOR PORT LOCK PROCEDURES

To calculate a close approximation of port system volume, you will need to determine the length of catheter used for each individual patient. (For future reference, it will be helpful to record this information on the patient's chart.) For the SmartPort® or SmartPort Plastic device, use formula and tables below:

$$\text{Port System Volume} = \text{Catheter Cut Length in cm} \times \frac{\text{catheter volume}}{\text{cm}} + \text{reservoir volume}$$

| SmartPort® and SmartPort Plastic device | |
|---|-------------------------------|
| Reservoir volume | 0.7 mL |
| Catheter | Catheter volume cm |
| 5F (1.8 mm) SL | 0.011 mL/cm |
| 6F (2.1 mm) SL | 0.013 mL/cm |
| 8F (2.6 mm) SL | 0.021 mL/cm |

If the port catheter length is not known, or generalized flushing volumes are desired per institutional protocols, the following flushing volumes are recommended:

| Procedure | Volume (100 U/mL) |
|--|--|
| When port not in use | 10 mL 0.9% Sodium Chloride every 4 weeks or based upon institutional policies and procedures |
| After each infusion of medication or TPN | 10 mL 0.9% Sodium Chloride for Injection |
| After blood withdrawal | 20 mL 0.9% Sodium Chloride for Injection |
| After Power Injection of Contrast Media | 10 mL 0.9% Sodium Chloride for Injection |

| | |
|--------------------------|---|
| Locking Procedure | To lock with heparin, repeat with 5 mL 100 U/mL heparinized saline, or with volume calculated above after use or based upon institutional policies and procedure. Close clamp while injecting last 0.5 mL of lock solution. |
|--------------------------|---|

Equipment ⁶

- Non-coring needle
- 10 mL syringe filled with 0.9% Sodium Chloride for Injection
- 10 mL syringe filled with 5 mL heparinized saline (100 U/mL)

Procedure ⁶

- Explain procedure to patient and prepare injection site.
- Perform hand hygiene before and after all vascular access procedures.
- Attach a 10-mL syringe filled with 0.9% Sodium Chloride for Injection to non-coring needle.
- Prepare skin with antiseptic solution per institutional protocol.
- Aseptically locate and access port including use of sterile gloves and mask per the non-coring needle DFU.
- Assess port functionality by using a 10 ml syringe or a syringe specifically designed to generate lower injection pressure (i.e. 10 ml diameter syringe barrel), taking note of any resistance to aspiration or flushing. Refer to System Occlusion Management section if resistance is encountered.
- Add dressing and stabilization per institutional protocol.
- Assess the access site, port functionality, change dressing and replace non-coring needle per institutional protocols.
- After therapy completion, flush port per institutional protocol. To lock with heparin, repeat with 5 mL 100 U/mL heparinized saline, or with volume calculated above. Close clamp while injecting last 0.5 mL of lock solution after use or based upon institutional policies and procedure.

- When therapy is complete, remove non-coring needle per product DFU and cover site according to institutional protocol.

Precaution: Do not use alcohol to soak or decontaminate polyurethane catheters because alcohol is known to degrade the polyurethane material over time with repeated and prolonged exposure.

SYSTEM OCCLUSION MANAGEMENT

Lumen obstruction is usually evident by failure to aspirate or infuse through the lumen or inadequate flow and/or high resistance pressures during aspiration and/or infusion. The causes may include inadequate catheter tip position, catheter kink, or catheter/vessel thrombosis, or fibrin sheath formation. If difficulty is encountered administering fluids through the SmartPort® or SmartPort Plastic device, discontinue the procedure and consider the following causes before proceeding with fibrinolytic therapy:

- Check non-coring needle placement. The non-coring needle may not be fully inserted through the septum of the port. Re-access the port and advance the needle until it is in contact with the bottom of the port chamber.
- Since blockage can be due to system movement, instruct the patient to alter position (e.g., move arms up and down, sit or stand).
- Have the patient cough.
- Provided there is no resistance, flush the catheter with 0.9% Sodium Chloride for Injection to try to move the tip away from the vessel wall.
- If suspected that blockage is due to blood clot in catheter, follow hospital/institutional policies. Use a 10 mL or larger syringe for resolution of blockage.

Precaution: Never forcibly flush an obstructed lumen.

If the lumen develops a thrombus, first attempt to aspirate the clot with a syringe. If aspiration fails, the clinician may consider using a thrombus dissolving solution to dissolve the clot. Refer to manufacturer's instructions, indications and precautions. If above methods prove unsuccessful, x-ray should be taken, as additional procedures may be required.

DISCONTINUING USE

If it is determined that a SmartPort® or SmartPort Plastic device is no longer required for therapy, the clinician should consider explanation of the system. If the system is left in place, then periodic x-ray with the patient in an upright, arms-at-the-sides position is recommended to detect problems with the system, such as pinching of the catheter between the clavicle and first rib which could result in catheter fragmentation and subsequent embolization.

WARRANTY

AngioDynamics warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond AngioDynamics control directly affect the instrument and the results obtained from its use. AngioDynamics obligation under this warranty is limited to the repair or replacement of this instrument and AngioDynamics shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. AngioDynamics neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **AngioDynamics assumes no liability with respect to instruments reused, reprocessed, resterilized, modified or altered in any way and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

REFERENCES

1. Stevens, Barbara, et al. "A Randomized, Prospective Trial of Conventional Vascular Ports vs. the Vortex 'Clear-Flow' Reservoir Port in Adult Oncology Patients." *Journal of Vascular Access Devices*, vol. 5, no. 2, 2000, pp. 37–40.
2. Chopra, Vineet, et al. "The Michigan appropriateness guide for intravenous catheters (Magic): Results from a multispecialty panel using the RAND/UCLA appropriateness method." *Annals of Internal Medicine*, vol. 163, no. 6_Supplement, 15 Sept. 2015, <https://doi.org/10.7326/m15-0744>
3. Shaw, Colette M., et al. "ACR appropriateness criteria ® radiologic management of Central Venous Access." *Journal of the American College of Radiology*, vol. 14, no. 11, Nov. 2017, <https://doi.org/10.1016/j.jacr.2017.08.053>
4. AngioDynamics Post Market Clinical Follow-up data on file for port survival and catheter related blood stream infection.
5. Hinke, D.H.; Zandt-Stastny, D.A.; Goodman, L.R.' et al. Pinch-off Syndrome: A Complication of Implantable Subclavian Venous Access Devices. *Radiology* 177: 353-356, 1990.
6. Kirmse, Chelsea Backler, J. *Access Device Guidelines: Recommendations for Nursing Practice and Education (Fourth Edition)*. Available from: Oncology Nursing Society, (4th Edition). Oncology Nursing Society, 2023.

* AngioDynamics, the AngioDynamics logo, SmartPort®, the SmartPort® logo, SmartPort, the SmartPort logo, Smart Angle, Vortex, BioFlo and LifeGuard are trademarks and/or registered trademarks of AngioDynamics, Inc., an affiliate or a subsidiary. All other trademarks are property of their respective owners.

Symbols for US: In compliance with the requirements of 21CFR Part 801.15, a glossary of symbols which appear without accompanying text within the product labeling is provided below

| Symbol | Ref | Title of Symbol | Meaning of Symbol |
|---|--------|---|--|
|  | 5.1.1 | Manufacturer | Indicates the medical device manufacturer. ^a |
|  | 5.1.2 | Authorized representative in the European Community / European Union. | Indicates the authorized representative in the European Community / European Union. ^a |
|  | 5.1.3 | Date of manufacture | Indicates the date when the medical device was manufactured. ^a |
|  | 5.1.4 | Use-by date | Indicates the date after which the medical device is not to be used. ^a |
|  | 5.1.5 | Batch code | Indicates the manufacturer's batch code so that the batch or lot can be identified. ^a |
|  | 5.1.6 | Catalog number | Indicates the manufacturer's catalog number so that the medical device can be identified. ^a |
|  | 5.1.8 | Importer | Indicates the entity importing the medical device into the locale. ^a |
|  | 5.2.3 | Sterilized using ethylene oxide | Indicates the medical device has been sterilized using ethylene oxide. ^a |
|  | 5.2.6 | Do not resterilize | Indicates a medical device that is not to be resterilized. ^a |
|  | 5.2.8 | Do not use if package is damaged and consult instructions for use | Indicates that a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information. ^a |
|  | 5.2.12 | Double sterile barrier system | Indicates a double sterile barrier system. ^a |
|  | 5.3.2 | Keep away from sunlight | Indicates a medical device that needs protection from light sources. ^a |
|  | 5.3.4 | Keep dry | Indicates a medical device that needs to be protected from moisture. ^a |

| Symbol | Ref | Title of Symbol | Meaning of Symbol |
|---|--------|--|--|
| | 5.4.2 | Do not re-use | Indicates a medical device that is intended for one single use only, or for use on a single patient during a single procedure. ^a |
| | 5.4.3 | Consult Instructions for use or consult electronic instructions for use ifu.angiodynamic.com | Indicates the need for the user to consult the instruction for use. ^a |
| | 5.4.10 | Contains hazardous substances | Indicates a medical device that contains substances that can be carcinogenic, mutagenic, reprotoxic (CMR) or substances with endocrine disrupting properties. ^a Contains cobalt as a component of stainless steel at levels $\leq 0.4\%$. This device is not intended for use in the stomach. Exposure of the stainless steel to highly acidic fluids such as gastric fluid can result in leaching of the cobalt from the stainless steel. Cobalt is listed in EC 1272/2008 as a carcinogen class 1B and a reproductive toxin class 1B. |
| | 5.7.7 | Medical device | Indicates the items is a medical device. ^a |
| | 5.7.10 | Unique Device Identifier | Indicates a carrier that contains unique device identifier information. ^a |
| | NA | Rx only | Caution: (US) Federal law restricts this device to sale by or on the order of a licensed practitioner. ^b |
| | NA | Universal Product Number | A Universal Product Number (UPN) code represents the manufacturer's number for an item. |
| | NA | Quantity in package | To indicate that the adjacent number reflects the number of units contained in the package. |
| | NA | CE Mark | Manufacturer's declaration of conformity to the Medical Device Regulation EU 2017/745. ^c |
| | 3.1.11 | Magnetic Resonance Conditional | An item with demonstrated safety in the MR environment within defined conditions including conditions for the static magnetic field, the time-varying gradient magnetic fields and the radiofrequency fields. |
| | NA | Contrast Enhanced Computed Tomography | Contrast Enhanced Computed Tomography |
| | 1135 | Recyclable Package | Recyclable Package. ^e |
| <p>a. EN ISO 15223-1 - Medical Devices – Symbols to be used with medical device labels, labeling and information to be supplied. b. 21 CFR 801.109 - Code of Federal Regulations. c. EU 2017/745 Medical Device Regulations published 5 May 2017 d. ASTM F2503-20 - Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment e. EN ISO 14021 Environmental labels and declarations. Self-declared environmental claims (Type II environmental labeling)</p> | | | |

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1. Stevens, Barbara, et al. "A Randomized, Prospective Trial of Conventional Vascular Ports vs. the Vortex 'Clear-Flow' Reservoir Port in Adult Oncology Patients." *Journal of Vascular Access Devices*, vol. 5, no. 2, 2000, pp. 37–40.
2. Chopra, Vineet, et al. "The Michigan appropriateness guide for intravenous catheters (Magic): Results from a multispecialty panel using the RAND/UCLA appropriateness method." *Annals of Internal Medicine*, vol. 163, no. 6_Supplement, 15 Sept. 2015, <https://doi.org/10.7326/m15-0744>
3. Shaw, Colette M., et al. "ACR appropriateness criteria® radiologic management of Central Venous Access." *Journal of the American College of Radiology*, vol. 14, no. 11, Nov. 2017, <https://doi.org/10.1016/j.jacr.2017.08.053>
4. AngioDynamics Post Market Clinical Follow-up data on file for port survival and catheter related blood stream infection.
5. Hinke, D.H.; Zandt-Stastny, D.A.; Goodman, L.R.' et al. Pinch-off Syndrome: A Complication of Implantable Subclavian Venous Access Devices. *Radiology* 177: 353-356, 1990.
6. Kirmse, Chelsea Backler, J. Access Device Guidelines: Recommendations for Nursing Practice and Education (Fourth Edition). Disponible de: Oncology Nursing Society, (4th Edition). Oncology Nursing Society, 2023.

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1. Stevens, Barbara, et al. "A Randomized, Prospective Trial of Conventional Vascular Ports vs. the Vortex 'Clear-Flow' Reservoir Port in Adult Oncology Patients." *Journal of Vascular Access Devices*, vol. 5, no. 2, 2000, pp. 37–40.
2. Chopra, Vineet, et al. "The Michigan appropriateness guide for intravenous catheters (Magic): Results from a multispecialty panel using the RAND/UCLA appropriateness method." *Annals of Internal Medicine*, vol. 163, no. 6_Supplement, 15 Sept. 2015, <https://doi.org/10.7326/m15-0744>
3. Shaw, Colette M., et al. "ACR appropriateness criteria® radiologic management of Central Venous Access." *Journal of the American College of Radiology*, vol. 14, no. 11, Nov. 2017, <https://doi.org/10.1016/j.jacr.2017.08.053>
4. AngioDynamics Post Market Clinical Follow-up data on file for port survival and catheter related blood stream infection.
5. Hinke, D.H.; Zandt-Stastny, D.A.; Goodman, L.R.' et al. Pinch-off Syndrome: A Complication of Implantable Subclavian Venous Access Devices. *Radiology* 177: 353-356, 1990.
6. Kirmse, Chelsea Backler, J. Access Device Guidelines: Recommendations for Nursing Practice and Education (Fourth Edition). Διαθέσιμο από: Oncology Nursing Society, (4th Edition). Oncology Nursing Society, 2023.

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1. Stevens, Barbara, et al., „A Randomized, Prospective Trial of Conventional Vascular Ports vs. the Vortex ‘Clear-Flow’ Reservoir Port in Adult Oncology Patients”. Journal of Vascular Access Devices, vol. 5, no. 2, 2000, pp. 37–40.
2. Chopra, Vineet, et al., „The Michigan appropriateness guide for intravenous catheters (Magic): Results from a multispecialty panel using the RAND/UCLA appropriateness method”. Annals of Internal Medicine, vol. 163, no. 6_Supplement, 15 Sept. 2015, <https://doi.org/10.7326/m15-0744>
3. Shaw, Colette M., et al., „ACR appropriateness criteria® radiologic management of Central Venous Access”. Journal of the American College of Radiology, vol. 14, no. 11, Nov. 2017, <https://doi.org/10.1016/j.jacr.2017.08.053>
4. AngioDynamics Post Market Clinical Follow-up data on file for port survival and catheter related blood stream infection.
5. Hinke, D.H.; Zandt-Stastny, D.A.; Goodman, L.R.' et al. Pinch-off Syndrome: A Complication of Implantable Subclavian Venous Access Devices. Radiology 177: 353–356, 1990.
6. Kirmse, Chelsea Backler, J. Access Device Guidelines: Recommendations for Nursing Practice and Education (Fourth Edition). Dokument dostępny w: Oncology Nursing Society, (4th Edition). Oncology Nursing Society, 2023.

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BIBLIOGRAFIA

1. Stevens, Barbara, et al. "A Randomized, Prospective Trial of Conventional Vascular Ports vs. the Vortex 'Clear-Flow' Reservoir Port in Adult Oncology Patients." *Journal of Vascular Access Devices*, vol. 5, no. 2, 2000, pp. 37–40.
2. Chopra, Vineet, et al. "The Michigan appropriateness guide for intravenous catheters (Magic): Results from a multispecialty panel using the RAND/UCLA appropriateness method." *Annals of Internal Medicine*, vol. 163, no. 6_Supplement, 15 Sept. 2015, <https://doi.org/10.7326/m15-0744>
3. Shaw, Colette M., et al. "ACR appropriateness criteria® radiologic management of Central Venous Access." *Journal of the American College of Radiology*, vol. 14, no. 11, Nov. 2017, <https://doi.org/10.1016/j.jacr.2017.08.053>
4. AngioDynamics Post Market Clinical Follow-up data on file for port survival and catheter related blood stream infection.
5. Hinke, D.H.; Zandt-Stastny, D.A.; Goodman, L.R.' et al. Pinch-off Syndrome: A Complication of Implantable Subclavian Venous Access Devices. *Radiology* 177: 353-356, 1990.
6. Kirmse, Chelsea Backler, J. Access Device Guidelines: Recommendations for Nursing Practice and Education (Fourth Edition). Disponibile da: Oncology Nursing Society, (4th Edition). Oncology Nursing Society, 2023.

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