VenaCure EVLT

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16900393-01 Rev. A

2022-10

NOTE: Read this manual and the Laser Operator's manual thoroughly prior to using the VenaCure EVLT* procedure kit. Observe all warnings, precautions and cautions noted throughout these documents. Failure to do so may result in patient complications.

INTENDED USE/ INDICATIONS FOR USE:

The AngioDynamics® VenaCure EVLT NeverTouch® Procedure Kits are indicated for endovascular coagulation of the Great Saphenous Vein (GSV) in patients with superficial vein reflux, for the treatment of varicose veins and varicosities associated with superficial reflux of the Great Saphenous Vein (GSV), and for the treatment of incompetence and reflux of superficial veins of the lower extremity.

This product should be used only with lasers cleared for use in the treatment of varicose veins, varicosities with superficial reflux of the GSV, and in the treatment of incompetent refluxing veins in the superficial venous system in the lower limbs.

INTENDED USER PROFILE

The VenaCure EVLT NeverTouch Procedure Kit is intended to be used by medical professionals, including but not limited to vascular surgeons, interventional cardiologists, and interventional radiologists. The use of the endovenous laser ablation devices must be restricted to medical staff with experience in the specialties suitable for the Indications for Use. It is recommended that all staff in contact with the device is informed about all applicable safety rules and standards ahead of use.

INTENDED PATIENT POPULATION

The VenaCure EVLT NeverTouch Procedure Kit is intended for use in patient populations that may include a broad range of ages, weights, races, nationalities, general health, and medical conditions (in accordance with all Warnings, Precautions, and Contraindications).

The VenaCure EVLT NeverTouch Procedure Kit is intended to treat patients with varicose veins.

CLINICAL BENEFITS

The expected clinical outcome of endovenous laser ablation (inclusive of the VenaCure EVLT NeverTouch Procedure Kit) is successful vein closure with minimal to no complications, and improvement in pain and clinical symptoms such as oedema, skin changes, inflammation and ulcers (as measured by the Venous Clinical Severity Score (VCSS)).

EXCLUSIONS:

- Patients with thrombus in the vein segment to be treated
- Patients with an aneurysmal section in the vein segment to be treated
- Patients with peripheral artery disease as determined by an Ankle-Brachial Index <0.9
- Patients with an inability to ambulate.
- Patients with deep vein thrombosis (DVT).
- Patients who are pregnant or breast feeding.
- Patients in general poor health.
- Other contraindications may be raised by the individual physician at the time of treatment.

POTENTIAL COMPLICATIONS:

The potential complications include but are not limited to the following: vessel perforation, thrombosis, pulmonary embolism, phlebitis, hematoma, infection, skin pigmentation alteration, neovascularization, paresthesia due to thermal damage of adjacent sensory nerves, anesthetic tumescence, non-target irradiation, hemorrhage, necrosis, air embolism, scarring, skin burns and pain.

PRECAUTION:

- Prior to and during use, avoid damaging the fiber by striking, stressing, or excessive bending. Do not coil the fiber tighter than a diameter of 20cm. Clinical safety and effectiveness data is not available for other fiber tip designs and diameters.
- Prior to and during use, avoid bending the introducer sheath and dilator as this can cause kinks and damage.
- All persons in the treatment room must be wearing protective glasses with the proper rating for the wavelength being used.

WARNINGS:

- Treatment of a vein located close to the skin surface may result in skin burn.
- Paresthesia may occur from thermal damage to adjacent sensory nerves.
- Tissue not targeted for treatment must be protected from injury by direct and reflected laser energy. Use appropriate eye and protective wear for both patient and operating personnel.
- Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your sales representative. Inspect prior to use to verify that no damage has occurred during shipping.
- Storage instructions: Keep away from sunlight; Keep dry; Upper limit of temperature: 27 degrees Celsius
- Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness, or death of the patient.
- Reprocessing may compromise the integrity of the device and/or lead to device failure, which, in turn, may result in patient injury, illness or death.
- VenaCure EVLT NeverTouch Procedure Kit is to be treated as contaminated biomedical waste subsequent to use. Additionally, the needle provided within the kit is a sharp device. The used or unused devices, as well as the sharps, should be disposed of in accordance with hospital, administrative and/or local government policy for such items.
- Uncontaminated device packaging should be recycled if applicable, or disposed of as common waste in accordance with hospital, administrative and/or local government policy for such items.
- The fiber, guidewire and needle contain cobalt. Cobalt is classified as CMR 1B and is present in a concentration above 0.1% weight by weight.

INSTRUMENTS REQUIRED:

- Laser
- Safety eyewear

NEVERTOUCH KIT CONTENTS:

- NeverTouch Gold-Tipped Fiber with SMA 905 connector and Sheath-Lok® Fitting
- 4 French Trè-Sheath® Introducer and Dilator
- 0.018" and/or 0.035" guidewire
- Needle
- 5 French, 10cm sheath/dilator - (optional)

INDIVIDUALLY PACKAGED Trè-Sheath CONTENTS:

- 4 French Trè-Sheath Introducer and Dilator

PROCEDURE:

Carefully read all instructions and observe all Warnings prior to performing the procedure. Patient complications may occur if this is not done. This procedure can be performed in the physicians' office or as an outpatient treatment under local anesthesia and should be performed by a qualified physician who has received training in these techniques.

1. Carry out a physical examination.
2. Review Ultrasound vein diameter and depth measurements.
3. Mark the skin with a surgical pen to show the length of vein to be treated.
4. Set the laser up according to the provided "Setting up the Laser Procedure" operators manual.

5. Place the patient in reverse Trendelenburg position so as to engorge and distend the targeted treatment vein.
6. Using standard Seldinger technique gain access to the desired vein with the needle and corresponding guidewire.
7. If applicable, use the micro-access introducer to exchange for a larger diameter wire.
8. Advance the 4Fr Trè-Sheath introducer, over the wire, to the treatment location. If treating the Great Saphenous Vein, position sheath tip so that it is 1-2cm below the Sapheno-Femoral Junction. Verify sheath positioning using ultrasound guidance.
9. Remove the dilator and guidewire.
10. Flush the sheath with normal saline.
11. Properly attach the sterile, single use disposable NeverTouch laser fiber by removing the protective cap and completely screwing the fiber fitting to lasers' SMA 905 connector until tight and secure.
12. Insert the NeverTouch laser fiber into the sheath until the Fiber Stop Assembly comes in contact with the end of the Sheath Hub. (See below)

13. Remove Fiber Stop Assembly. Pull the sheath back and lock it to the Sheath-Lok Fitting.
14. Confirm location of the fiber using ultrasound guidance.
15. The NeverTouch fiber end should be 2 cm below the Sapheno-Femoral junction when treating the GSV.
16. Move the Sliding Sheath Gauge until flush with the insertion site, to mark the final position of the sheath.
17. Administer local tumescent anesthesia under ultrasound guidance as to ensure delivery of just enough anesthesia to achieve thermal protection.
18. Set the laser to continuous mode and adjust to desired power.
(810nm, 980nm Lasers: 10-14 Watts)
(1470nm Lasers: 5-7 Watts)
19. Put on laser safety glasses appropriate for the applicable wavelength.
20. Place the laser in the Enable mode.
21. Check position of fiber end and sliding sheath gauge; adjust position if necessary.
22. Activate the laser by depressing the foot pedal while withdrawing the fiber and sheath together, at a rate that adequately delivers desired laser energy per centimeter.
(810nm, 980nm Lasers: 50-80 Joules per cm)
(1470nm Lasers: 30-50 Joules per cm)
Do not apply direct external hand pressure or force over the fiber tip during energy activation.
23. Cease operating the laser by removing foot from foot pedal when the fiber end is 2 – 3 cm from the access site as indicated by markers on the distal tip of the Trè-Sheath Introducer. The procedure is considered complete when the desired length of the targeted treatment vein is treated.

POST PROCEDURE:

Perform a second ultrasound to confirm that the vein is closed and that there is no more blood flow.

Any serious incident which has occurred with the use of this device should be reported to AngioDynamics or your in-country representative.

Any serious incident which has occurred with the use of this device should be reported to AngioDynamics at complaints@angiodynamics.com and to the National Competent Authority. Refer to the following web address for contact information of the EU Competent Authority. https://ec.europa.eu/health/sites/health/files/md_sector/docs/md_vigilance_contact_points.pdf
The instructions for use are available electronically at www.angiodynamics.com/ifu-dfu-portal/


WARRANTY

AngioDynamics, Inc. warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other representations and warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond AngioDynamics, Inc.'s control directly affect the instrument and the results obtained from its use. AngioDynamics, Inc.'s obligation under this warranty is limited to the repair or replacement of this instrument and AngioDynamics, Inc. shall not be liable for any incidental or consequential losses, damages or expenses directly or indirectly arising from the use of this instrument. AngioDynamics, Inc. neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. AngioDynamics, Inc. assumes no liability with respect to instruments reused, reprocessed or resterilized, modified or altered in any way, and makes no warranties or representations, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.

* AngioDynamics, the AngioDynamics logo, NeverTouch, VenaCure EVLT, the VenaCure EVLT logo, Trè-Sheath, and Sheath-Lok are trademarks and/or registered trademarks of AngioDynamics Inc., an affiliate or a subsidiary.

In compliance with the requirements of 21 CFR Part 801.15, below is a glossary of symbols which appear without accompanying text.

Symbol	Ref	Title of Symbol	Meaning of Symbol
	5.1.1	Manufacturer	Indicates the medical device manufacturer. ^a
	5.1.2	Authorized representative in the European Community/ European Union	Indicates the authorized representative in the European Community/European Union. ^a
	5.1.3	Date of manufacture	Indicates the date when the medical device was manufactured. ^a

Symbol	Ref	Title of Symbol	Meaning of Symbol
	5.1.4	Use-by date	Indicates the date after which the medical device is not to be used. ^a
LOT	5.1.5	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified. ^a
REF	5.1.6	Catalog number	Indicates the manufacturer's catalog number so that the medical device can be identified. ^a
	5.1.8	Importer	Indicates the entity importing the medical device into the locale. ^a
STERILE EO	5.2.3	Sterilized using ethylene oxide	Indicates the medical device has been sterilized using ethylene oxide. ^a
	5.2.6	Do not re-sterilize	Indicates a medical device that is not to be re-sterilized. ^a
	5.2.8	Do not use if package is damaged and consult instructions for use	Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information. ^a
	5.2.11	Single sterile barrier system	Indicates a single sterile barrier system. ^a
	5.3.2	Keep away from sunlight	Indicates a medical device that needs protection from light sources. ^a
	5.3.4	Keep dry	Indicates a medical device that needs to be protected from moisture. ^a
	5.3.6	Upper limit of temperature	Indicates the upper limit of temperature to which the medical device can be safely exposed. ^a
	5.4.2	Do not re-use	Indicates a medical device that is intended for one single use only, or for use on a single patient during a single procedure. ^a
	5.4.3	Consult instructions for use or consult electronic instructions for use ifu.angiodynamics.com	Indicates the need for the user to consult the instruction for use. ^a
 Cobalt	5.4.10	Contains hazardous substances	Indicates a medical device that contains substances that can be carcinogenic, mutagenic, reprotoxic (CMR) or substances with endocrine disrupting properties. ^a Contains cobalt as a component of stainless steel. This device is not intended for use in the stomach. Exposure of the stainless steel to highly acidic fluids such as gastric fluid can result in leaching of the cobalt from the stainless steel. Cobalt is listed in EC 1272/2008 as a carcinogen class 1B and a reproductive toxin class 1B.
MD	5.7.7	Medical device	Indicates the items is a medical device. ^a
UDI	5.7.10	Unique Device Identifier	Indicates a carrier that contains unique device identifier information. ^a
Rx ONLY	NA	Rx only	Caution: (US) Federal law restricts this device to sale by or on the order of a licensed practitioner. ^b
UPN	NA	Universal Product Number	A Universal Product Number (UPN) code represents the manufacturer's number for an item.
	NA	Quantity in package	To indicate that the adjacent number reflects the number of units contained in the package.
CE 2797	NA	CE Mark	Manufacturer's declaration of conformity to the Medical Device Regulation EU 2017/745. ^c
	1135	Recyclable Package	Recyclable Package. ^d

Symbol	Ref	Title of Symbol	Meaning of Symbol
a.		EN ISO 15223-1 - Medical Devices – Symbols to be used with medical device labels, labeling and information to be supplied.	
b.		21 CFR 801.109 - Code of Federal Regulations.	
c.		EU 2017/745 Medical Device Regulations published 5 May 2017	
d.		EN ISO 14021 Environmental labels and declarations. Self-declared environmental claims (Type II environmental labeling)	

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AngioDynamics UK Ltd.