# ANGIOGRAPHIC CATHETER

## C € 2797

### Recommended Instructions for Use

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16900444-01A  
2022-08

## INDICATIONS

ANGIODYNAMICS® Angiographic Catheters are for use where angiographic diagnosis is indicated.

These instructions for use are available electronically at [ifu.angiodynamics.com](http://ifu.angiodynamics.com)

### Intended Use

Angiographic catheters are single-use invasive medical devices which are intended to be used for patients in need of angiographic diagnosis. The catheter is inserted into the vasculature system and guided to the target anatomical location.

The primary function of AngioDynamics angiographic catheters is to deliver contrast media to a specified anatomical location for subsequent diagnosis or intervention.

Any serious incident which has occurred with the use of this device should be reported to AngioDynamics or your in-country representative .

### Notice for European Union Only:

- For a copy of the Summary of Safety and Clinical Performance for this device, please review Eudamed ([ec.europa.eu/tools/eudamed](http://ec.europa.eu/tools/eudamed)) referencing UDI-DI # 505168400552 or contact AngioDynamics Customer Service at +1 800-772-6446
- Any serious incident which has occurred with the use of this device should be reported to AngioDynamics at [complaints@angiodynamics.com](mailto:complaints@angiodynamics.com) and to the National Competent Authority. Refer to the following web address for contact information for the Competent Authorities. [https://ec.europa.eu/health/sites/health/files/md\\_sector/docs/md\\_vigilance\\_contact\\_points.pdf](https://ec.europa.eu/health/sites/health/files/md_sector/docs/md_vigilance_contact_points.pdf)
- The physician should discuss the risks associated with the device with the patient.

### Product Description

The AngioDynamics Angiographic Catheters are provided in a variety of tip shapes and endhole configurations to suit the physician's needs depending on the region of the vasculature where the procedure is being performed. Unique features within the product families include:

- Braided or non-braided co-extruded nylon shaft with Soft, radiopaque tips for easy visualisation under fluoroscopy
- Sizing catheter with consistently placed radiopaque marker bands across the entire sizing pattern (within +/- 1mm accuracy) (Accu-Vu only).
- Hydrophilic coating to the distal portion of the catheter reducing catheter surface friction (Mariner only)

### Patient target groups

Angiographic catheters are intended to be used for patients that require subsequent diagnosis, intervention, or evaluation of blood vessels throughout the body.

### Expected Clinical Benefits

Angiographic catheters are single-use invasive medical devices which are intended to be used for patients in need of angiographic diagnosis. The catheter is inserted into the vessel and guided through various arteries to the target organ. The primary function is to deliver contrast media to a specified vessel for subsequent diagnosis or intervention.

The Intended Clinical Benefit of AngioDynamics Angiographic Catheters is the delivery of contrast media for the purpose of angiographic imaging and examination. The contrast medium, a radiopaque dye, allows visualization of the target organ through x-rays, subsequently allowing for the identification of diseases of blood vessels such as:

- Aneurysm: a weakening of the arterial wall
- Arteriosclerosis: hardening of the arteries, loss of elasticity of the arterial wall
- Atherosclerosis: fatty substance on the inner wall of an artery reducing its lumen and restricting blood flow
- Embolus: a traveling blood clot or foreign material
- Stenosis: constriction or narrowing of the vessel
- Occlusion: blockage of the vessel
- Tumor: abnormal mass

Specific to the hydrophilic coating on the *Mariner*® Hydrophilic-Coated Angiographic Catheter, the hydrophilic coating applied to the distal portion of the outer shaft makes the catheter feel slippery to the touch when it is wetted and is intended to enhance catheter trackability by reducing surface friction. This feature is intended to facilitate smooth insertion and withdrawal of the catheter within the vasculature to prevent vessel injury.

### POTENTIAL ADVERSE EFFECTS

The following adverse reactions have been reported and are associated with the use of angiographic catheters:

- |                        |                       |                         |
|------------------------|-----------------------|-------------------------|
| • Thrombus formation   | • Plaque dislodgment  | • Myocardial infarction |
| • Emboli               | • Hematoma            | • Stroke                |
| • Arterial wall damage | • Cardiac arrhythmias | • Death                 |

### WARNINGS

1. ANGIODYNAMICS Angiographic Catheters should be used only by physicians with a thorough understanding of angiography and percutaneous interventional procedures. This includes Interventional Radiologists, Vascular Surgeons, Interventional Cardiologists and other Physicians with a thorough understanding of angiography and percutaneous interventional procedures.
2. Do not insert catheters directly through synthetic vascular grafts. Insert through a sheath introducer.
3. Never advance or retract an Angiographic Catheter or guidewire against resistance. This may result in damage to the vessel, the product, or both.
4. Do not attempt to hand straighten the tip of any curved tipped catheters which are furnished with a tip straightener. Hand straightening may result in damage to the product. Tip straighteners are furnished on all catheters intended to be straightened with the aid of a tip straightener.

5. Reshaping of the catheter tip is not recommended. Physical damage to the catheter material can result when exposed to heat.
6. Do not leave curved tipped catheters straightened over guidewires for extended periods of time. This will result in failure of the catheter tip to re-form to its intended shape.
7. When using ANGIODYNAMICS HALO Catheters, form or straighten the tip in the distal thoracic aorta using a guidewire when advancing or withdrawing the catheter.
8. Avoid injection of contrast medium into the intercostal and lumbar arteries.
9. Always use a guidewire to remove the catheter from the vasculature. Failure to do so may result in damage to the vessel, puncture site, product, or all three.
10. When retracting catheters, great care must be taken to avoid exerting excessive pressure at the groin entry site. Excessive pressure may result in damage to the vessel, the catheter, or both.
11. The maximum pressure limit of AngioDynamics Angiographic catheters intended for flush angiography is stated on the catheter package. When using a pressure injector, do not exceed the stated maximum pressure.
12. AngioDynamics Angiographic Catheters intended for selective angiography do not have a maximum pressure limit stated on the catheter package. Typical flow rates of up to 10cc per second are stated with pressure generated to achieve these typical flow rates.
13. Flow rates have been established by using a contrast medium as stated on the catheter label at 37° centigrade for a 2 second injection. Since a number of variables (such as media viscosity and temperature) can produce different flow rates, this information should only be used as a guide.
14. Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your AngioDynamics, Inc. representative. Inspect prior to use to verify that no damage has occurred in shipping.
15. For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
16. AngioDynamics Angiographic Catheters are to be treated as contaminated biomedical waste subsequent to use. Used or unused devices should be disposed of in accordance with hospital, administrative and/or local government policy for such items.
17. Uncontaminated device packaging should be recycled if applicable or disposed of as common waste in accordance with hospital, administrative and/or local government policy for such items.

#### INSTRUCTIONS FOR USE

- A. Introduce the Angiographic Catheter into the vasculature using appropriate vascular access technique for the patient.
- B. Precautions to reduce clotting should be taken when any catheter is used. Use of systemic heparinization should be considered especially in small vessels.
- C. Utilize a guidewire for advancement, repositioning or removal of the catheter.
  - a. Refer to the AngioDynamics catheter labeling for the recommended guidewire diameter to be used with the catheter for optimal performance. Guidewire use has been associated with a greater incidence of thrombus formation. Optimal guidewire size and judicious use are recommended.
  - b. Refer to guidewire manufacturer's instructions for use for proper prep and use.
  - c. If resistance is encountered when removing the guidewire from the catheter, simultaneously remove the catheter and guidewire as a single unit under fluoroscopy to prevent potential vessel and product damage.
- D. For catheters fitted with a tip straightener, utilize the device to straighten out the catheter tip. Do not straighten tip using other means.
  - a. To straighten the tip, hold the distal end of the tip straightener between the thumb and forefinger and then pull the straightener to the distal end of the catheter.
  - b. Once the tip is straight, place the catheter tip over the positioned guidewire and introduce the catheter into the vasculature.
  - c. To remove the tip straightener from the catheter, gently slide the straightener to the proximal end of the catheter (towards the strain relief). Continue to pull gently until the straightener opens and slides off the catheter.
- E. While in the vascular system, the catheter should always be filled with either contrast medium or saline.
- F. Prior to injections using a pressure injector, always confirm the position of the distal tip of the catheter using fluoroscopic guidance. Check for free backflow of blood and inject a small amount of contrast medium by hand if necessary, to confirm the catheter tip's position.
  - a. To check for catheter's tip position, prior to injection, check for the backflow of blood and inject a small amount of contrast medium by hand
- G. Prior to use and injections, refer to the catheter packaging for guidance on:
  - a. Selective catheters: the generated pressures at common flow rates
  - b. Flush catheters: pressure limit of 1050psi (7239 kpa) and the maximum flow rate (ml/sec). Do not exceed the maximum pressure.
- H. If resistance is encountered during catheter manipulation, determine the cause under fluoroscopy and take the necessary remedial action.
- I. For AngioDynamics Accu-Vu sizing catheters, the marker band measurement accuracy is  $\pm 1$ mm.

## WARRANTY

AngioDynamics, Inc. warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other representations and warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond AngioDynamics, Inc.'s control directly affect the instrument and the results obtained from its use. AngioDynamics, Inc.'s obligation under this warranty is limited to the repair or replacement of this instrument and AngioDynamics, Inc. shall not be liable for any incidental or consequential losses, damages or expenses directly or indirectly arising from the use of this instrument. AngioDynamics, Inc. neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **AngioDynamics, Inc. assumes no liability with respect to instruments reused, reprocessed or resterilized, modified or altered in any way, and makes no warranties or representations, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

\* AngioDynamics, the AngioDynamics logo, Accu-Vu, Soft-Vu and Mariner are trademarks and/or registered trademarks of AngioDynamics, Inc.

In compliance with the requirements of 21 CFR Part 801.15, below is a glossary of symbols which appear without accompanying text.

Symbol	Ref	Title of Symbol	Meaning of Symbol
	5.1.1	Manufacturer	Indicates the medical device manufacturer. <sup>a</sup>
	5.1.2	Authorized representative in the European Community	Indicates the authorized representative in the European Community. <sup>a</sup>
	5.1.3	Date of manufacture	Indicates the date when the medical device was manufactured. <sup>a</sup>
	5.1.4	Use-by date	Indicates the date after which the medical device is not to be used. <sup>a</sup>
	5.1.5	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified. <sup>a</sup>
	5.1.6	Catalog number	Indicates the manufacturer's catalog number so that the medical device can be identified. <sup>a</sup>
	5.1.8	Importer	Indicates the entity importing the medical device into the locale. <sup>a</sup>
	5.2.3	Sterilized using ethylene oxide	Indicates the medical device has been sterilized using ethylene oxide. <sup>a</sup>
	5.2.6	Do not re-sterilize	Indicates a medical device that is not to be re-sterilized. <sup>a</sup>
	5.2.8	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened. <sup>a</sup>
	5.2.11	Single sterile barrier system	Indicates a single sterile barrier system. <sup>a</sup>
	5.3.2	Keep away from sunlight	Indicates a medical device that needs protection from light sources. <sup>a</sup>
	5.3.4	Keep dry	Indicates a medical device that needs to be protected from moisture. <sup>a</sup>
	5.3.6	Upper limit of temperature	Indicates the upper limit of temperature to which the medical device can be safely exposed. <sup>a</sup>
	5.4.2	Do not re-use	Indicates a medical device that is intended for one single use only, or for use on a single patient during a single procedure. <sup>a</sup>
	5.4.3	Consult instructions for use <a href="http://ifu.angiodynamics.com">ifu.angiodynamics.com</a>	Indicates the need for the user to consult the instruction for use. <sup>a</sup>
	5.7.7	Medical device	Indicates the items is a medical device. <sup>a</sup>
	5.7.10	Unique Device Identifier	Indicates a carrier that contains unique device identifier information. <sup>a</sup>

<b>Rx ONLY</b>	NA	Rx only	Caution: (US) Federal law restricts this device to sale by or on the order of a licensed practitioner. <sup>b</sup>
<b>UPN</b>	NA	Universal Product Number	A Universal Product Number (UPN) code represents the manufacturer's number for an item.
	NA	Quantity in package	To indicate that the adjacent number reflects the number of units contained in the package.
<b>CE 2797</b>	NA	CE Mark	Manufacturer's declaration of conformity to the Medical Device Regulation EU 2017/745. <sup>c</sup>
	NA	Maximum Guidewire OD	Maximum Guidewire Outer Diameter
	1135	Recyclable Package	Recyclable Package. <sup>d</sup>
<p>a. EN ISO 15223-1 - Medical Devices – Symbols to be used with medical device labels, labeling and information to be supplied.</p> <p>b. 21 CFR 801.109 - Code of Federal Regulations.</p> <p>c. EU 2017/745 Medical Device Regulations published 5 May 2017</p> <p>d. EN ISO 14021 Environmental labels and declarations. Self-declared environmental claims (Type II environmental labeling)</p>			

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AngioDynamics UK Ltd.