angiodynamics

NanoKnife

EN	SINGLE ELECTRODE ACTIVATION PROBE INSTRUCTIONS FOR USE	2
DE	EINZELEKTRODEN-AKTIVIERUNGSSONDE GEBRAUCHSANWEISUNG	6
ES	SONDA DE ACTIVACIÓN DE ELECTRODO ÚNICO INSTRUCCIONES DE USO	11
FR	SONDE D'ACTIVATION À ÉLECTRODE SIMPLE MODE D'EMPLOI.....	16
IT	SONDA DI ATTIVAZIONE A ELETTRODO SINGOLO ISTRUZIONI PER L'USO	20
NL	GEBRUIKSAANWIJZING AKTIVERINGSSONDE MET ÉÉN ELEKTRODE	24
DA	AKTIVERINGSSONDE MED EN ELEKTRODE BRUGSANVISNING	28
SV	SOND MED EN ELEKTROD OCH AKTIVERINGSSOND BRUKSANVISNING.....	32
NO	AKTIVERINGSSONDE MED ENKELTELEKTRODE BRUKSANVISNING.....	36
PT	SONDA DE ATIVAÇÃO DE ELÉTRODO INDIVIDUAL INSTRUÇÕES DE UTILIZAÇÃO.....	40

AngioDynamics, Inc.
603 Queensbury Avenue
Queensbury, NY 12804
USA
USA Customer Service 800-772-6446

CE 2797

EC REP

AngioDynamics Netherlands BV
Haaksbergweg 75
1101 BR, Amsterdam
The Netherlands

**UK
CA
0086**

UK Responsible Person:
AngioDynamics UK Ltd.
c/o Kidd Rapinet
29 Harbour
Exchange Square
London E14 9GE
United Kingdom

WARNINGS

Do NOT use a device with damaged insulation.

Do NOT attach anything to the device unless it is supplied by AngioDynamics and indicated for use with this device. Attachments may damage the insulation and contribute to patient injury.

MR Unsafe: This device should not be used in conjunction with MRI image guidance as the device has not been tested for MRI compatibility.

CAUTION

Rx ONLY Federal (USA) law restricts this device to sale by or use under the order of a physician.

This device is sterilized by ethylene oxide and intended for single patient use. Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient. Store the NanoKnife Single Electrode Probe in a cool, dry place. Inspect prior to use to verify no damage has occurred during shipping. Do not use if sterile barrier is received open or damaged. Do not re-sterilize. Reprocessing may compromise the integrity of the device and / or lead to device failure. Do not use if the labeling is incomplete or illegible.

- NanoKnife probes are sharp devices. Used and unused devices should be disposed in accordance with hospital, administrative, and/or local government policy for such devices
- NanoKnife probe spacers, after use, are contaminated biomedical waste. The used device should be disposed in accordance with hospital, administrative, and/or local government policy for such waste.
- The probe is intended for use with an *AngioDynamics *NanoKnife generator only. The physician must read the NanoKnife System User Manual thoroughly before operating the unit.
- Avoid having the NanoKnife Single Electrode Probes touch when delivering pulses. Ensure that both NanoKnife Single Electrode Probe electrodes are fully embedded within the targeted tissue prior to initiating pulse delivery.
- Do not bend the NanoKnife Single Electrode Probe as this may damage the insulation.
- Do not attach anything to the device unless it is supplied by AngioDynamics and indicated for use with this device. Attachments may damage the insulation and contribute to patient injury. Placing a Single Electrode Probe without continual image guidance may increase the risk of unintended mechanical perforation, damage to critical anatomical structures, and/or unintended tissue destruction. Imaging equipment with high definition capabilities is recommended for procedures where a target area involves or is adjacent to critical structures. Having inadequate imaging equipment may result in unintended mechanical perforation, damage to critical anatomical structures, and/or hemorrhage.
- Refer to the NanoKnife Generator User Manual for a complete list of warnings and precautions.

POTENTIAL ADVERSE EFFECTS

Adverse effects that may be associated with the use of the NanoKnife System include, but are not limited to:

- Arrhythmia
 - Atrial fibrillation or flutter
 - Bigeminy
 - Bradycardia
 - Heart block or atrioventricular block
 - Paroxysmal supraventricular tachycardia
 - Tachycardia
 - > Reflex tachycardia
 - > Ventricular tachycardia
 - Ventricular fibrillation
- Damage to critical anatomical structure (nerve, vessel, duct)
- Fistula formation
- Hematoma
- Hemorrhage
- Hemothorax
- Infection
- Muscle contraction
- Pneumothorax
- Reflex hypertension
- Unintended mechanical perforation
- Vagal stimulation, asystole
- Venous thrombosis

INTENDED USE

The Single Electrode Activation Probes are intended for use with the NanoKnife System.

PRODUCT DESCRIPTION

The Single Electrode Activation Probes have 19 gauge trocar needles and are available in 15 cm and 25 cm lengths. Including the insulation, diameter is 18 gauge.


Figure 1: Probe configurations

The probes are comprised of:

1. Active electrode, length adjustable in 0.5 cm increments from 0 – 4 cm via the thumbslide.
2. Insulation sleeve.
3. 19 gauge needle with depth markers and echogenic needle tip.
4. 10 foot (3 m) connection cable.
5. Sterile labels for identifying the generator port number the probe will be connected to.


Figure 2: Probe components

The Single Electrode Activation Probe may be attached to any output port of the generator. Depending on the size of the soft tissue area to be ablated, a maximum of six probes can be used at a time. Probes may be re-positioned after each ablation to cover a larger area, as guided by the generator software.

INSTRUCTIONS FOR USE

1. Refer to the NanoKnife System User Manual and become familiar with operation of the generator.
2. Turn on the generator.
3. Inspect all of the devices and packaging for damage prior to use. Do NOT use any devices that are damaged or if the sterile barrier is breached. Unused probes from an opened procedure pack should be discarded after the procedure.
4. Using sterile technique, open the package and carefully remove the device. Remove and dispose of the protective shipping sheath covering the needle.
5. Sterile labels are included with the Single Electrode Activation Probe that may be used to identify the connector and corresponding handles that will be attached to each port of the generator during the procedure.


Figure 3: Labels for connector and handle

6. Attach the cable connectors to the front panel of the generator. **Note:** If excessive force is required to connect the cable to the generator, check to see if the pin in the cable connector is bent.
7. Program the generator with the parameters required to achieve the desired ablation zone. Refer to the generator manual for specific programming instructions.
 - a. Typical results (based on in vivo porcine data) for two probes spaced 1 cm apart with the active electrode length set at 2 cm.

Generator Settings:

Probe +	Probe -	Voltage	Pulse Length	N. Pulses
1	2	2000 V	100 μ s	90


Figure 4: Typical ablation results based on in vivo porcine data.

- b. Typical results (based on in vivo porcine data) for two probes spaced 1.5 cm apart with the active electrode length set at 2 cm.

Generator Settings:

Probe +	Probe -	Voltage	Pulse Length	N. Pulses
1	2	2500 V	100 μ s	90


Figure 5: Typical ablation results based on in vivo porcine data.

8. Adjust the exposure of the electrode using the thumbslide, pressing down on the proximal raised end to unlock. Exposure may be adjusted from 0 to 4.0 cm in 0.5 cm increments. The accuracy of the exposure adjustment is +/- 0.15 cm.


Figure 6: Electrode exposure adjustment

9. Verify entry point and trajectory of each Single Electrode Probe prior to insertion using imaging equipment. Using imaging guidance, the probes should be placed in the targeted tissue area as parallel as possible. Placing a Single Electrode Probe without continual image guidance may increase the risk of unintended mechanical perforation, damage to critical anatomical structures, and/or unintended tissue destruction. Use one or more NanoKnife Probe Spacers to help place the Single Electrode Probes parallel to each other and at a set distance apart.
10. Place each Single Electrode Probe, carefully and systematically, using imaging equipment, under continual image guidance, while avoiding tissue obstacles and critical structures. The probes should be placed such that they are as parallel as possible to one another. Avoid placing the Single Electrode Probes such that they converge at the tip. Converging probe placement may lead to high current during pulse delivery.
11. Markings on the needle can assist in placement of the device. The markings extend up the needle in 1.0 cm increments.

Figure 7: Depth markings on needle

12. Verify the position of the probes with imaging. Needles should be verified to be as close to parallel as possible and consistent with the probe placement plan.
13. Verify the length of the exposed electrode indicated by the thumb slide on the handle (see Figure 6).
14. Once all desired probes are placed, measure and record all the inter-electrode distances using imaging equipment measurement tools.
15. Refer to generator manual for pulse delivery instructions.

PROBE REMOVAL AND DISPOSAL

16. After the ablation is complete before removing the probes, review the graphs on the Pulse Generation screen to ensure that the pulses were delivered.
17. Use the thumb slide to change the Probe Exposure Setting to 0 cm for each probe, thereby covering the exposed electrode region.
18. Remove each probe from the patient and place onto the sterile prep area.
19. Follow medical judgement to achieve hemostasis at puncture site after probe removal.
20. Disconnect all NanoKnife probe cable connectors from the NanoKnife Generator
21. NanoKnife probes are sharp devices. Used and unused devices should be disposed in accordance with hospital, administrative, and/or local government policy for such devices
22. NanoKnife probe spacers, after use, are contaminated biomedical waste. The used device should be disposed in accordance with hospital, administrative, and/or local government policy for such waste.

23. Untampered device packaging should be recycled if applicable, or disposed of as common waste in accordance with hospital, administrative and/or local government policy for such items.

* AngioDynamics, the AngioDynamics logo, NanoKnife and the NanoKnife logo are trademarks and/or registered trademarks of AngioDynamics, Inc., an affiliate or a subsidiary.

In compliance with the requirements of 21 CFR Part 801.15, below is a glossary of symbols which appear without accompanying text.

Symbol	Ref	Title of Symbol	Meaning of Symbol
	5.1.1	Manufacturer	Indicates the medical device manufacturer. ^a
	5.1.2	Authorized representative in the European Community/ European Union	Indicates the authorized representative in the European Community/European Union. ^a
	5.1.3	Date of manufacture	Indicates the date when the medical device was manufactured. ^a
	5.1.4	Use-by date	Indicates the date after which the medical device is not to be used. ^a
	5.1.5	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified. ^a
	5.1.6	Catalog number	Indicates the manufacturer's catalog number so that the medical device can be identified. ^a
	5.1.8	Importer	Indicates the entity importing the medical device into the locale. ^a
	5.2.3	Sterilized using ethylene oxide	Indicates the medical device has been sterilized using ethylene oxide. ^a
	5.2.6	Do not resterilize	Indicates a medical device that is not to be resterilized. ^a
	5.2.8	Do not use if package is damaged and consult instructions for use	Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information. ^a
	5.2.11	Single sterile barrier system	Indicates a single sterile barrier system. ^a
	5.3.2	Keep away from sunlight	Indicates a medical device that needs protection from light sources. ^a
	5.3.4	Keep dry	Indicates a medical device that needs to be protected from moisture. ^a
	5.3.6	Upper limit of temperature	Indicates the upper limit of temperature to which the medical device can be safely exposed. ^a
	5.4.2	Do not re-use	Indicates a medical device that is intended for one single use only, or for use on a single patient during a single procedure. ^a
	5.4.3	Consult instructions for use or consult electronic instructions for use ifu.angiodynamics.com	Indicates the need for the user to consult the instruction for use. ^a
	5.7.7	Medical device	Indicates the items is a medical device. ^a
	5.7.10	Unique Device Identifier	Indicates a carrier that contains unique device identifier information. ^a
	NA	Rx only	Caution: (US) Federal law restricts this device to sale by or on the order of a licensed practitioner. ^a
	NA	Universal Product Number	A Universal Product Number (UPN) code represents the manufacturer's number for an item.
	NA	Quantity in package	To indicate that the adjacent number reflects the number of units contained in the package.
	NA	CE Mark	Manufacturer's declaration of conformity to the Medical Device Regulation EU 2017/745. ^c
	NA	Magnetic Resonance (MR) unsafe	Keep away from magnetic resonance imaging (MRI) equipment. ^a
	1135	Recyclable Package	Recyclable Package. ^a

a. EN ISO 15223-1 - Medical Devices – Symbols to be used with medical device labels, labeling and information to be supplied.

b. 21 CFR 801.109 - Code of Federal Regulations.

c. EU 2017/745 Medical Device Regulations published 5 May 2017

d. ASTM F2503-13 - Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.

e. EN ISO 14021 Environmental labels and declarations. Self-declared environmental claims (Type II environmental labeling)

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